

Article code	2847001
Article name acc. DoC	Setback valve
Class acc. Directive 93/42/EEC	Ila

Article code	Article name	Material
C10745	Check valve, small	●reducing part: PC (Polycarbonate) ●connection part: PC (Polycarbonate) ●valve: LDPE(low density polyethylene)
C11515	Solvent	●MEK (Methylethylketone)
C11553	Solvent	●CHN (Cyclohexanone)
C22108	Cap for check valve small	●PP (Polypropylene)
C22502	N-Heptan for Spectroscopy	●Heptane
C30959	Tube 5.1 x 7.1 mm	●PVC (Polyvinylchloride) with plasticizer DEHP, Colourant blue
C30999	Silicone Tube 5 x 1,5 mm	●Silikon-Elastomer
C31525	Printing ink, blue	●Mixture (colourant and thinners)
C31526	Thinner for printing ink	●1-Methoxy-2-Propanol (2-Metholymethylethoxy) propanol
Packaging		
	White paper for IFU printing	●100% wood pulp
	Self-adhesive tape, white	●PP (Polypropylene) ●adhesive: natural rubber
	Natron Paper 1000 x 750 mm	●Color brown
	Label sticker 160 x 115 mm	●Paper
	Stretch film	●LLDPE (linear low density polyethylene)
	Stretch film foil plate	●LDPE (low density polyethylene)
	Cardboard 586 x 251 x 165 mm	●Corrugated cardboard, color brown ●Printing FK Logo: Blue- Pantone 300U
	Protective Layer 200 x 370 mm	●cardboard
	Compound foil 30/70 µm	●PA (Polyamide)/PE (Polyethylene)
	Automat reel without print	●Paper
Sterilization		
	Indicator tape for EO sterilization	●Ethylene Oxide (EO) sterilization
	Sterilization gas	●Ethylene oxide