

Article code	2867001
Article name acc. DoC	Freka Lipoflow
Class acc. Directive 93/42/EEC	Is

Article code	Article name	Material
C10463	Roberts clamp, white	●ABS (Acrylonitrile Butadiene Styrene)
C10942	Tube 3.0 x 4.1 mm	●Polyether Polyurethane
C10942A	Polyurethane Tube Ø3,00 x Ø4,1 mm	●Polyether Polyurethane
C11503	Solvent	●THF (Tetrahydrofuran), also 1,4-Epoxy butane
C21656	Withdrawal spike, vented, with protective cap	● Protective cap: LDPE (low density polyethylene) ● Spike: ABS (Acrylonitrile Butadiene Styrene) ● Filter holder: HDPE (Polyethylene) ● Filter membrane: Acrylic copolymer on non-woven nylon
C21657	Disposable Hypodermic Needle 16Gx25mm	●Needle tubing: stainless steel ●Needle hub: PP (Polypropylene) ●Protector: PP (Polypropylene) ●Lubricant: PDMS (Polydimethylsiloxane, Silicone oil) ●Glue: Epoxy-resin glue
C30179	Male connector LL	●MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene)
Packaging		
	Carton 400 x 220 x 380 mm	●Corrugated cardboard, Flute B
	Tubular foil for packaging	●Film 500m x 350 mm x 40 µm ●Mixture of different Polyolefines PE (Polyethylene), LDPE (Polyethylene, low density), transparent
	Compound foil 70 µm	●PA (Polyamide) / adhesive layer / PE (Polyethylene)
	PE banderole 80 x 30 x 0.10 mm	●LDPE (low density polyethylene)
	Label sticker, size: 104 x 148 mm	●Paper for laserprinting
	White paper for printing, IFU	●100% wood pulp
	Adhesive tape	●BOPP (biaxially oriented polypropylene) film and rubber adhesive
	Paper Freka Lipoflow 100 x 180 mm	●Paper
	Stretch film	●LLDPE (linear low density polyethylene)
	Stretch film foil plate	●LDPE (low density polyethylene)
Sterilization		
	Indicator tape for EO sterilization	●Process indicator tape for Ethylene Oxide (EO) sterilization
	Sterilization gas	●100% Ethylene oxide