

Article code	2883101
Article name acc. DoC	Transfer Set 1 D
Class acc. MDR 2017/745	Ila

Article code	Article name	Material
C11515	Solvent	●MEK (Methylethylketone)
C11553	Solvent	●CHN (Cyclohexanone)
C21915	Tube adaptor piece	●TPE resp. SEBS (Styrene-Ethylene-Butylene-Styrene)
C60951	PVC Tube 5,10 x 7,10	●Plasticized Polyvinyl Chloride, plasticiser DEHT
J1010481		
C10176	Drip chamber upper part with spike	●ABS (Acrylonitrile-Butadiene-Styrene) ●Colourant white
C11515	Solvent	●MEK (Methylethylketone)
C11553	Solvent	●CHN (Cyclohexanone)
C22112	Protective cap	●LDPE (low density polyethylene)
C22116	Drip chamber lower part	●PVC (Polyvinylchloride)
C22559	Glue	●ECA (Ethyl-2-Cyanoacrylate) ; Hydrochinon
J1010851		
C22113	Roller clamp housing	●PC (Polycarbonate) ●Colourant blue, base resin PC (Polycarbonate)
C22114	Roller clamp wheel	●ABS (Acrylonitrile Butadiene Styrene) ●Colourant white
J1040014		
C22105	Urology connector	●PVC (Polyvinylchloride), colourant white, plasticizer DINP
C22106	Protective cap fur urology adaptor small	●LDPE (low density polyethylene) + Masterbatch (Blue)
Packaging		
	Self-adhesive tape, white	●PP (Polypropylene) ●adhesive: natural rubber
	Natron Paper 1000 x 750 mm	●Color brown
	Label sticker 160 x 115 mm	●Paper
	Stretch film	●LLDPE (linear low density polyethylene)
	Stretch film foil plate	●LDPE (low density polyethylene)
	Cardboard 586 x 251 x 165 mm	●Corrugated cardboard, color brown ●Printing FK Logo: Blue- Pantone 300U
	Protective Layer 200 x 370 mm	●cardboard
	Compound foil 30/70 µm	●PA (Polyamide)/PE (Polyethylene)
	Printing ink, blue	●Mixture (colourant and thinners)
	Thinner for printing ink	●1-Methoxy-2-Propanol (2-Metholymethylethoxy) propanol

Article code	Article name	Material
	Automat reel without print	●Paper
Sterilization		
	Indicator tape for EO sterilization	●Process indicator tape for Ethylene Oxide (EO) sterilization
	Sterilization gas	●Ethylene oxide