

Article code	2886452
Article name acc. DoC	Infudrop Air-Nitro G, PVC-free, light-protected
Class acc. Directive 93/42/EEC	Is

Article code	Article name	Material
C10842	infusion flat filter 0.2 µm, positively charged	●Housing: polyacryl ●Filter: Supor® membrane, pore size 0.2 µm, hydrophilic polyethersulfone ●Vent membrane: 0.02 µm supported by PTFE (Polytetrafluoroethylene)
C10845	filter (tape down clip = back of C10842)	●modified polyacryl
C11503	Solvent	●THF (Tetrahydrofuran), also 1,4-Epoxy butane
C11515	Solvent	●MEK (Methylethylketone)
C30947	Tube 3 x 4.1 mm	●PUR (Polyurethane) + yellow UV stabilized Masterbatch PUR yellow
J1010112		
C10178	Drip chamber upper part 21 E	●PS (polystyrene) ●Colourant white, 50 % TiO₂ in HIPS (high impact polystyrene)
C10221	moulded ring white for drip chamber 21	●PS (Polystyrene) ●Colourant, carrier material PS (Polystyrene), contains 50% TiO₂
C10266	Drip chamber lower part 21 V92, yellow	●SBS (Styrene Butadiene Copolymer) ●masterbatch with UV-stabilizer, 5% yellow, Carrier: PS (Polystyrene)
C10345	air ventilation flap	●1-propene homopolymer: PP (Polypropylene) ●1-propene copolymer (heterophasic copolymer)
C10539	Protective cap	●PP (Polypropylene)
C10836	ventilation filter 3µm	●Membrane Type: Hydrophobic ●Material Type: ePTFE (expanded Polytetrafluorethylene) Membrane ●Membrane Support Type: Non-woven Polyester
J1010725		
C10533	protective cap	●PP (polypropylene) ●LDPE (low density polyethylene) with inorganic color pigments, white
C10762	Back check valve	●Valve membrane: Silicone ●lower part (Luer Lock male): SAN (Styrene acrylonitrile copolymer), coulorant white ●upper Part (tube fitment): MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene), colourant white ●Solvent: THF (Tetrahydrofuran)
J1010805		
C10404	roller clamp housing, white	●PS (polystyrene) ●PS (polystyrene) with silicone

Article code	Article name	Material
C10412	Roller clamp wheel	●PS (polystyrene) ●PS (Polystyrene) with silicone slip-masterbatch
Packaging		
	Carton box	●Cardboard
	Compound foil 70 µm	●PA (Polyamide) / adhesive layer / PE (Polyethylene)
	Adhesive tape paper	●gummed paper tapes - Polyeste filament yarn for strengthening
	PE-Band 90 x 30 x 0.05 mm	●HDPE (high density polyethylene)
	Label sticker, size: 104 x 148 mm	●Paper for laserprinting
	White paper for IFU printing	●100% wood pulp
	Stretch film	●LLDPE (linear low density polyethylene)
	Stretch film foil plate	●LDPE (low density polyethylene)
	Paper Infudrop light protected 100 x 225 mm	●Paper
J1010112		
	PE bag	●LDPE (Low density polyethylene)
Sterilization		
	Indicator tape for EO sterilization	●Process indicator tape for Ethylene Oxide (EO) sterilization
	Sterilization gas, Poland	●Ethylene oxide