

Article code	2890701
Article name acc. DoC	Twin Transfer Set 9 D
Class acc. MDR 2017/745	Ila

Article code	Article name	Material
C11515	Solvent	●MEK (Methylethylketone)
C11553	Solvent	●CHN (Cyclohexanone)
C21913	Tube clamp	●POM (Polyoxymethylene) or Acetal-Copolymer, base resin, additives ●Resin ester polymer (carrier material) + pigment dispersion
C21915	Tube adaptor piece	●TPE resp. SEBS (Styrene-Ethylene-Butylene-Styrene)
C30723	Y-piece, parallel	●PVC (Polyvinylchloride) noDOP ●Colourant blue
C60946	PVC Tube 6,20 x 8,80	●Plasticized Polyvinyl Chloride, plasticiser DEHT
C60951	PVC Tube 5,10 x 7,10	●Plasticized Polyvinyl Chloride, plasticiser DEHT
J1040053	Piercing spike with grip plate	●N/A
J1010483		
C11515	Solvent	●MEK (Methylethylketone)
C11553	Solvent	●CHN (Cyclohexanone)
C22116	Drip chamber lower part	●PVC (Polyvinylchloride)
C22123	Drip chamber upper part	●PVC (Polyvinylchloride) shore 65 D with DEHT (Di-2-ethylhexyl-terephthalate) plasticizer
J1010851		
C22113	Roller clamp housing	●PC (Polycarbonate) ●Colourant blue, base resin PC (Polycarbonate)
C22114	Roller clamp wheel	●ABS (Acrylonitrile Butadiene Styrene) ●Colourant white
J1040014		
C22105	Urology connector	●PVC (Polyvinylchloride), colourant white, plasticizer DINP
C22106	Protective cap for urology adaptor small	●LDPE (low density polyethylene) + Masterbatch (Blue)
Packaging		
	Self-adhesive tape, white	●PP (Polypropylene) ●adhesive: natural rubber
	Natron Paper 1000 x 750 mm	●Color brown
	Label sticker 160 x 115 mm	●Paper
	Stretch film	●LLDPE (linear low density polyethylene)
	Stretch film foil plate	●LDPE (low density polyethylene)

Article code	Article name	Material
	Cardboard 586 x 251 x 165 mm	●Corrugated cardboard, color brown ●Printing FK Logo: Blue- Pantone 300U
	Protective Layer 200 x 370 mm	●cardboard
	Compound foil 30/70 µm	●PA (Polyamide)/PE (Polyethylene)
	Printing ink, blue	●Mixture (colourant and thinners)
	Thinner for printing ink	●1-Methoxy-2-Propanol (2-Metholymethylethoxy) propanol
	Automat reel without print	●Paper
Sterilization		
	Indicator tape for EO sterilization	●Process indicator tape for Ethylene Oxide (EO) sterilization
	Sterilization gas, Poland	●Ethylene oxide