

Article code	2892100
Article name acc. DoC	Ambix activ set ambulatory SF
Class acc. Directive 93/42/EEC	Ila

Article code	Article name	Material
C10564	Protective cap LL pos. with 1.2 µm filter	●Housing: PP (Polypropylene) + Masterbatch ●Filter: Haemophobic/hydrophobic membrane; 1,2 µm PVDF (Polyvinylidene difluoride) Fortex Membrane
C10747	Back check valve	●Housing upper part: MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene), white ●Housing lower part: SAN (Styrene-Acrylonitrile) ● Colourant: white ● Valve: Silicone rubber
C10898	Filter Speed Flow 1,2 µm	●Housing: MBS (Methylacrylate / butadiene/styrene) ●fluid filter: 1.2 µm hydrophilic PES (Polyethersulfone) ●air filter: 0.03 µm hydrophobic PTFE (Polytetrafluoroethylene)
C11503	Solvent	●THF (Tetrahydrofuran), also 1,4-Epoxy butane
C11553	Solvent	●CHN (Cyclohexanone)
C21650	Non-vented spike	●Protective Cap: PE (Polyethylene) ●Spike: ABS (Acrylonitrile Butadiene Styrene), white
C30949	Tube 2.5 x 4.22 mm	●PVC (Polyvinylchloride)
J1010953		
C10464	Clamp	● PC (Polycarbonate)
C11503	Solvent	●THF (Tetrahydrofuran), also 1,4-Epoxy butane
C11515	Solvent	●MEK (Methylethylketone)
C30631	Connection piece	● MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene) ● (2%) carrier material ABS (Acrylonitrile Butadiene Styrene) with anorganic and organic colour pigments (blue)
C30632	Pump connector upstream	● MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene) ● (2%) carrier material PS (Polystyrene) with anorganic and organic pigments and dyes (light green)
C30650	Antikink	●PP (Polypropylene) ●(2%) Carrier material PP with organinc and inorganic colourants
C30957	Tube 3.0 x 5.0 x 167 ± 5 mm	●PVC (Polyvinylchloride) noDOP, with plasticizer TEHTM (Tri-(2-ethylhexyl) trimellitate)
Packaging		
	Folding box 400 x 220 x 95 mm	●Corrugated cardboard
	Compound foil 70 µm	●PA (Polyamide) / adhesive layer / PE (Polyethylene)

Article code	Article name	Material
	PE-Band 90 x 30 x 0.05 mm	●HDPE (high density polyethylene)
	White paper for IFU printing	●100% wood pulp
	Roll Label 105 x 148 mm	●Paper (blank white) and Adhesive
	Stretch film	●LLDPE (linear low density polyethylene)
	Stretch film foil plate	●LDPE (low density polyethylene)
	Paper "Ambix activ Set"	●Paper
Sterilization		
	Indicator tape for EO sterilization	●Ethylene Oxide (EO) sterilization
	Sterilization gas	●Ethylene oxide