

Article code	2894201
Article name acc. DoC	Twin Transfer Set 1 D
Class acc. MDR 2017/745	Ila

Article code	Article name	Material
C10744	Check valve, large	<ul style="list-style-type: none"> ●Housing: HDPE (high density Polyethylene) ●Valve plates: EPDM (Ethylene Propylene Diene Rubber)
C11515	Solvent	●MEK (Methylethylketone)
C11553	Solvent	●CHN (Cyclohexanone)
C21879	Connecting funnel	<ul style="list-style-type: none"> ●TPU (thermoplastic Polyurethane) ●Preparation of a colour agent in a universal carrier ●Lubricant (Amid wax) ●TPU-(thermoplastic Polyurethane)-based anti-blocking agent
C21915	Tube adaptor piece	●TPE resp. SEBS (Styrene-Ethylene-Butylene-Styrene)
C30723	Y-piece, parallel	<ul style="list-style-type: none"> ●PVC (Polyvinylchloride) noDOP ●Colourant blue
C30963	Tube 7.98 x 2.1 mm	●PVC (Polyvinylchloride) with plasticizer
C60946	PVC Tube 6,20 x 8,80	●Plasticized Polyvinyl Chloride, plasticiser DEHT
C60951	PVC Tube 5,10 x 7,10	●Plasticized Polyvinyl Chloride, plasticiser DEHT
J1040053	Piercing spike with grip plate	●N/A
J1010851		
C22113	Roller clamp housing	<ul style="list-style-type: none"> ●PC (Polycarbonate) ●Colourant blue, base resin PC (Polycarbonate)
C22114	Roller clamp wheel	<ul style="list-style-type: none"> ●ABS (Acrylonitrile Butadiene Styrene) ●Colourant white
Packaging		
	Self-adhesive tape, white	<ul style="list-style-type: none"> ●PP (Polypropylene) ●adhesive: natural rubber
	Natron Paper 1000 x 750 mm	●Color brown
	Label sticker 160 x 115 mm	●Paper
	Stretch film	●LLDPE (linear low density polyethylene)
	Stretch film foil plate	●LDPE (low density polyethylene)
	Cardboard 586 x 251 x 165 mm	<ul style="list-style-type: none"> ●Corrugated cardboard, color brown ●Printing FK Logo: Blue- Pantone 300U
	Protective Layer 200 x 370 mm	●cardboard
	Compound foil 30/70 µm	●PA (Polyamide)/PE (Polyethylene)
	Printing ink, blue	●Mixture (colourant and thinners)
	Thinner for printing ink	●1-Methoxy-2-Propanol (2-Metholymethylethoxy) propanol

Article code	Article name	Material
	Automat reel without print	●Paper
Sterilization		
	Indicator tape for EO sterilization	●Process indicator tape for Ethylene Oxide (EO) sterilization
	Sterilization gas, Poland	●Ethylene oxide