

Article code	2909201
Article name acc. DoC	Freka® Connect 20ml ENFit/ProNeo (30x1)
Class acc. Directive 93/42/EEC	IIa

Article code	Article name	Material
B-0075-4	Infusionsfilter RoweFil 0,2 µm	●Infusionsfilter (welded)
E-0068	Protective cap	●PE (Polyethylene)
K-0021	Gas indicator tape	●N/A
K-0057	Adaptor LL male	●MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene)
K-0087	Tube 2.5 x 4.1 x 100 ± 5 mm	●PVC (Polyvinylchloride) with plasticizer DEHP
K-0179	Protective cap for LL female	●ABS (Acrylonitrile Butadiene Styrene), white
K-0320	Positively loaded Nylon-Membrane 0,2	●N/A
K-0337	Solvent	●THF (Tetrahydrofuran), CAS-no. 109-99-9
K-1213	PTFE-Membran 0,03 µm	●N/A
Packaging		
	carton 290 x 190 x 235 mm	●Corrugated cardboard
	Foil bag 300 x 500 mm	●PE (Polyethylene)
	Paper 105 x 134 mm	●over foil
	Compound foil 70 µm	●PA (Polyamide) / PE (Polyethylene)
	Compound foil 70 µm	●PA (Polyamide) / PE (Polyethylene)
	Label sticker 74 x 105 mm	●Paper
	Label sticker 105 x 148 mm	●carton product label
	instruction for use	●Paper
Sterilization		
	Sterilization	ETO (Ethylene Oxide)