

Article code	2909301
Article name acc. DoC	Lipifil-Air, 24 h PVC-frei
Class acc. Directive 93/42/EEC	IIa

Article code	Article name	Material
B-0172-3	Filter (housing with yellow arrows)	●Casing: MABS (K-0370) ●Capillary pore filter: PET (Polyethyleneterephthalate) 1.2 µm (K-0358) ●Membran: PTFE (K-1213)
E-0068	Protective cap	●PE (Polyethylene)
E-0088-1	Adaptor LL male	●MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene)
K-0021	Gas indicator tape	●N/A
K-0179	Protective cap for LL female	●ABS (Acrylonitrile Butadiene Styrene), white
K-0287	Tube 3.0 x 4.1 x 100 ± 5 mm	●PUR (Polyurethane)
K-0337	Solvent	●THF (Tetrahydrofuran), CAS-no. 109-99-9
Packaging		
	carton 290 x 190 x 235 mm	●Corrugated cardboard
	Foil bag 300 x 500 mm	●PE (Polyethylene)
	Paper 105 x 134 mm	●over foil
	Compound foil 70 µm	●PA (Polyamide) / PE (Polyethylene)
	Compound foil 70 µm	●PA (Polyamide) / PE (Polyethylene)
	instruction for use	●Paper
	Label sticker 105 x 148 mm	●carton product label
Sterilization		
	Sterilization	●Ethylene oxide