

Article code	2909502
Article name acc. DoC	Infufil Air 96h, 0.2 µm, positive, 10 cm2
Class acc. MDR 2017/745	Ila

Article code	Article name	Material
C10542	Protective cap	●HDPE (high density polyethylene)
C10607	LL neg. cone with tube fitment	●PVC (Polyvinylchloride) without plasticizer
C10867	Speedflow 0.2 micron positive charged membrane	●Filter media: Hydrophilic PES (Polyethersulfone) membrane 0.2 µm ●Vent: Hydrophobic PTFE (Polytetrafluoroethylene) 0.03 µm ●Housing: Clear Modified Acrylic
C11503	Solvent	●THF (Tetrahydrofuran), also 1,4-Epoxy butane
C30936	Tube 3.0 x 4.1 mm	●PVC (Polyvinylchloride) soft, noDOP, 34% TEHTM (Tri-(2-ethylhexyl) trimellitate)
J1010725		
C10533	protective cap	●PP (polypropylene) ●LDPE (low density polyethylene) with inorganic color pigments, white
C10762	Back check valve	●Valve membrane: Silicone ●lower part (Luer Lock male): SAN (Styrene acrylonitrile copolymer), colourant white ●upper Part (tube fitment): MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene), colourant white ●Solvent: THF (Tetrahydrofuran)
Packaging		
	Carton 400x220x380 mm	●3-layer cardboard
	Compound foil 70 µm	●PA (Polyamide) / adhesive layer / PE (Polyethylene)
	Adhesive tape paper	●gummed paper tapes - Polyeste filament yarn for strengthening
	White paper for IFU printing	●100% wood pulp
	Roll Label 105 x 148 mm	●Paper (blank white) and Adhesive
	PE bag with bottomseal, transparent	●LDPE (low density polyethylene) ●Film thickness: 45 µm
	Stretch film	●LLDPE (linear low density polyethylene)
	Stretch film foil plate	●LDPE (low density polyethylene)
	PE (Polyethylene) film broadband	●PE (Polyethylene), white
	Peel paper	●Paper

Article code	Article name	Material
Sterilization		
	Indicator tape for EO sterilization	●Process indicator tape for Ethylene Oxide (EO) sterilization
	Sterilization gas	●Ethylene oxide