

Article code	2909602
Article name acc. DoC	Infufil Air 96h, 0.2 µm, 10cm ²
Class acc. MDR 2017/745	Ila

Article code	Article name	Material
C10542	Protective cap	●HDPE (high density polyethylene)
C10607	LL neg. cone with tube fitment	●PVC (Polyvinylchloride) without plasticizer
C10899	Speedflow adult filter 0,2µm	●Filter media: Hydrophilic PES (Polyethersulfone) membrane ●Vent: Hydrophobic PTFE (Polytetrafluoroethylene) 0.03 µm ●Housing: Clear Modified Acrylic
C11503	Solvent	●THF (Tetrahydrofuran), also 1,4-Epoxy butane
C30936	Tube 3.0 x 4.1 mm	●PVC (Polyvinylchloride) soft, noDOP, 34% TEHTM (Tri-(2-ethylhexyl) trimellitate)
J1010725		
C10533	protective cap	●PP (polypropylene) ●LDPE (low density polyethylene) with inorganic color pigments, white
C10762	Back check valve	●Valve membrane: Silicone ●lower part (Luer Lock male): SAN (Styrene acrylonitrile copolymer), coulorant white ●upper Part (tube fitment): MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene), colourant white ●Solvent: THF (Tetrahydrofuran)
Packaging		
	Carton 400 x 220 x 190 mm	●Cardboard
	Compound foil 70 µm	●PA (Polyamide) / adhesive layer / PE (Polyethylene)
	Adhesive tape paper	●gummed paper tapes Polyeste filament yarn for strengthening
	Label sticker, size: 104 x 148 mm	●Paper for laserprinting
	White paper for IFU printing	●100% wood pulp
	Roll Label 105 x 148 mm	●Paper (blank white) and Adhesive
	PE bag with bottomseal, transparent	●LDPE (low density polyethylene) ●Film thickness: 45 µm
	Stretch film	●LLDPE (linear low density polyethylene)
	Stretch film foil plate	●LDPE (low density polyethylene)

	PE (Polyethylene) film broadband	●PE (Polyethylene), white
	Peel paper	●Paper
	Peel paper	●Paper
Sterilization		
	Sterilization gas, Nanchang	●Ethylene oxide (90%) + Carbon dioxide (10%)
	Indicator tape for EO sterilization	●Process indicator tape for Ethylene Oxide (EO) sterilization
	Sterilization gas, Poland	●100% Ethylene oxide, liquefied grade one (GAZ 4-134) ●alternative: > 99.9% Ethylene oxide