

Article code	2909802
Article name acc. DoC	Infufil 5 µm, 4.6 cm ²
Class acc. Directive 93/42/EEC	Is

Article code	Article name	Material
C10833	Filter	●Housing: ABS ●Filter: Polyester
Packaging		
	Carton 400 x 220 x 380 mm	●Corrugated cardboard, Flute B
	Tubular foil for packaging	●Film 500m x 350 mm x 40 µm ●Mixture of different Polyolefines PE (Polyethylene), LDPE (Polyethylene, low density), transparent
	Compound foil 70 µm	●PA (Polyamide) / adhesive layer / PE (Polyethylene)
	Adhesive tape paper	●gummed paper tapes - Polyeste filament yarn for strengthening
	Label sticker, size: 104 x 148 mm	●Paper for laserprinting
	White paper for IFU printing	●100% wood pulp
	Stretch film	●LLDPE (linear low density polyethylene)
	Stretch film foil plate	●LDPE (low density polyethylene)
	Paper Streudr Fresenius	●Paper
Sterilization		
	Indicator tape for EO sterilization	●Ethylene Oxide (EO) sterilization
	Sterilization gas	●Ethylene oxide