

Material List
 Finished Product: 5036741



Internal code	Designation	Polymer
Material 224	Protection Cap	●PP (Polypropylene)
Material 46	Cannula straight, diameter 0.9mm, stitch length 37 mm	●Stainless steel
Material 54	Hub with LL	●PC (Polycarbonate) ●colourant yellow
Material 55	Glue for cannula	●Epoxy resin

Packaging

Material 56	Blister 102 x 27.85 x 14 mm	●PA (Polyamide)/PE (Polyethylene) foil 90 µm
Material 57	Lid 102 x 27.85 mm	●Bacteria-proof paper 80g/sqm
Material 58	Folding box 110 x 50 x 50 mm	●hard paper, white
Material 59	Label 67 x 38 mm	●Paper
Material 60	Carton 400 x 215 x 340 mm	●Corrugated cardboard, brown
Material 61	Sterile indicator	●N/A
Material 62	Tape	●PP (Polypropylene), transparent
Material 63	Stretch foil	●PE (Polyethylene), transparent
Material 64	EURO-Pallet 1200 x 800	●Wood