

Material List
 Finished Product: 5036771



Internal code	Designation	Polymer
Material 225	Tube 1x1,6	●N/A
Material 51	Cannula bent, diameter 0.9mm, stitch length 37 mm	●Stainless steel
Material 54	Hub with LL	●PC (Polycarbonate) ●colourant yellow
Material 55	Glue for cannula	●Epoxy resin

Packaging

Material 59	Label 67 x 38 mm	●Paper
Material 61	Sterile indicator	●N/A
Material 62	Tape	●PP (Polypropylene), transparent
Material 63	Stretch foil	●PE (Polyethylene), transparent
Material 64	EURO-Pallet 1200 x 800	●Wood
Material 65	Blister 67.5 x 65 x 9.5 mm	●PA (Polyamide)/PE (Polyethylene) foil 90 µm
Material 66	Lid 67.5 x 65 mm	●Bacteria-proof paper 80g/sqm
Material 67	Folding box 70 x70 x 70 mm	●hard paper, white
Material 68	Carton 380 x 380 x 375 mm	●Corrugated cardboard, brown