

Article code	7751011AP
Article name acc. DoC	APPLIX Gravity Set VarioLine
Class acc. Directive 93/42/EEC	Is

Article code	Article name	Material
C10467	Roller clamp	●housing: PP (Polypropylene) plus colourant blue ●wheel: HDPE (high density polyethylene)
C11198	Flat pouch 145 x 255 x 0.08 mm	●LDPE (Low density polyethylene)
C11503	Solvent	●THF (Tetrahydrofuran), also 1,4-Epoxy butane
C11553	Solvent	●CHN (Cyclohexanone)
C21815	Spike for spike port system	●ABS (Acrylonitrile-Butadiene-Styrene) ●ABS (Acrylonitrile-Butadiene-Styrene), Pigment white
C30059	Universal bottle adaptor	●LDPE (low density polyethylene) ●HDPE (high density polyethylene) ●Colourant white ●lubricant demoulding aid
C30558	ENLock universal male	●ABS (Acrylonitrile-Butadiene-Styrene) ●Pigment ABS purple
C41134	Drip chamber	●cap and lower part: SBC (Styrene Butadiene Copolymer) ●glue: THF (Tetrahydrofuran) ●adaptor for tube: PUR (Polyurethane)
J1010755		
C10330	Air ventilation plug	●ABS (Acrylonitrile Butadiene Styrene) ●Colourant white, ABS (Acrylonitrile Butadiene Styrene)-carrier
C10824	Bacteria filter	●glass microfiber filter
C11120	Lid valve	●Silicone
C30034	Universal bottle adapter	●LDPE (low density polyethylene) ●HDPE (high density polyethylene) ●LDPE (low density polyethylene) with inorganic color pigments, white ●lubricant demoulding aid
J1020016		
C10895	Sealing for Pouch connection	●PE (Polyethylene) foam
C21792	Retainer nut	●PP (Polypropylene) ●Colourant white
C21793	Spike for spike port system	●ABS (Acrylonitrile Butadiene Styrene) ●Colourant white
J1040004		

Article code	Article name	Material
C22189	Closure with adaptor cone	●POM (polyoxymethylene) Copolymer, granulate, base resin, additives ●POM (polyoxymethylene) masterbatch yellow
C22190	T-piece	● MABS (Methymethacrylate Acrylonitrile Butadiene Styrene)
J1040152		
C11503	Solvent	●THF (Tetrahydrofuran), also 1,4-Epoxy butane
C30606	Luer Lock adaptor female	●PVC (Polyvinylchloride)
C30953	Tube 3.5 x 5.0 mm	●PVC (Polyvinylchloride) with plasticizer TOTM
Packaging		
	Carton 1185 x 785	●Carton
	Tubular foil for packaging	●Film 500m x 350 mm x 40 µm ●Mixture of different Polyolefines PE (Polyethylene), LDPE (Polyethylene, low density), transparent
	Carton box	●Cardboard
	Carton box 391 x 288 x 209 mm	●Corrugated cardboard
	Adhesive tape paper	●gummed paper tapes - Polyeste filament yarn for strengthening
	PE-Band 90 x 30 x 0.05 mm	●HDPE (high density polyethylene)
	Label sticker, size: 104 x 148 mm	●Paper for laserprinting
	ENLock sticker	●80 g/m2 adhesive label, white, permanent
	Paper 100 x 225 mm	●Paper
	Cardboard plate 1170 x 770 mm	●Carton plate
	Cardboard plate 250 x 455 mm	●Carton plate
	PE bag with bottomseal, transparent	●LDPE (low density polyethylene) ●Film thickness: 45 µm
	Stretch film	●LLDPE (linear low density polyethylene)
	Stretch film foil plate	●LDPE (low density polyethylene)
	PE (Polyethylene) film broadband	●PE (Polyethylene), white
	Compound foil 90µm	●Multilayer flat film PA (Polyamide) /adhesive layer / PE (Polyethylene)
Sterilization		
	Indicator tape for EO sterilization	●Process indicator tape for Ethylene Oxide (EO) sterilization
	Sterilization gas	●Ethylene oxide