

Article code	7751102
Article name acc. DoC	Hydrobag
Class acc. Directive 93/42/EEC	Is

Component	Material	Plasticizer
Top filler	PVC DEHP-free	DOTP
Bag	PVC DEHP-free	DOTP
Tube	PVC DEHP-free	DOTP
Clamp	PP	DOTP
Ucomfor Connector	PVC DEHP-free	DOTP
Packaging		
Blister bag	Paper 63g PA/PE foil thickness: 0.16mm	-
IFU	Paper	-
Outer carton	Brown Corrugated paper 5layer	-
Barcode label	Adhesive Label	-
White paper card	Paper 350G	-
Sterilization		
	Indicator tape for EO sterilization	●Process indicator tape for Ethylene Oxide (EO) sterilization
	Sterilization gas	●100% Ethylene oxide