

Article code	7751531AP
Article name acc. DoC	Freka PEG Set gastric FR 20
Class acc. Directive 93/42/EEC	I Ib

Article code	Article name	Material
C21832	Spool and thread	● MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene), Colourant (blue) ● PA 6 / PA 6.6 (Polyamide)
C21837	Fixation screw	● PBT (Polybutyleneterephthalate) ● Resin ester polymer (carrier material) + pigment dispersion
C21839	Fixation ring Ø6,5 mm CH20	● TPE (Thermoplastic Elastomer)
C21843	Proximal cone	● TPU (thermoplastic Polyurethane), Polyester-based ● Resin ester polymer (carrier material) + pigment additive
C21845	Thread introducer device	● PP (Polypropylene) ● Resin ester polymer (carrier material) + pigment dispersion
C21857	Fixation aid PEG	● PBT (Polybutyleneterephthalate) ● Resin ester polymer (carrier material) + pigment dispersion
C21881	Connecting funnel with Luer	● TPU (thermoplastic Polyurethane), Polyether-based ● Resin ester polymer (carrier material) + pigment dispersion ● Lubricant (Amid wax)
C22606	Tray cover sheet	● coloured PETG (Polyethyleneterephthalate Glycol), blue tinted, mould release: silicone additive
C22609	PEG tray 200 x 155 mm	● coloured PETG (Polyethyleneterephthalate Glycol), blue tinted, mould release: silicone additive
C22663	Scalpel	● Blade: stainless steel ● Handle: HIPS (high impact polystyrene) ● Guard: LDPE (low density polyethylene)
C22691	Fixation loop, OD 0.4 ± 0.01mm	● Core: PA 6.6 (Polyamide 6.6) ● Coating: PA 6 (Polyamide 6)
C30992	Tube 5x6.6x350 ± 5mm	● PUR (Polyurethane) and PUR + 60% BaSO ₄ (one x-ray contrast strip)
C31531	Tube printing ink	● Ketone based colourant (cyclohexanone, butanone)
C31532	Tube printing ink thinner	● Ketone based additive (Butanone = Methyl ethyl ketone)
J1040029		
C11515	Solvent	● MEK (Methyl ethyl ketone)
C11553	Solvent	● CHN (Cyclohexanone)

Article code	Article name	Material
C21803	Blood stop plug	●MABS (Methymethacrylate Acrylonitrile Butadiene Styrene)
C22081	Valve, lower part	●PUR (Polyurethane) with masterbatch green and release additive
C22082	Cannula valve upper part	●PUR (Polyurethane) ●Colourant
C22084	Valve opening piece	●PP (Polypropylene)
C22085	Protective cap	●PE (Polyethylene)
C22558	Glue	●ECA (Ethyl-2-Cyanoacrylate); Hydrochinon
C22573	Lubricant	●Siloxane
C22588	Cannula	●needle holder: PC (polycarbonate based on BPA bisphenole A) ●needle: stainless steel ●glue Polyurethane
C22665	Capillary 8/80	●Fluorinated Ethylene Propylene (FEP); ●Stripes: Pantone 542 blue; ●X-ray contrast agent: 24% BaSO4 (pre-mixed in main material)
C22937S	Membrane	●LSR (Liquid silicone rubber)
C31757	PTFE (Polytetrafluoroethylene) Etching Agent	●Etching Agent for PTFE
J1040082		
C21835	PEG LL pos	●PBT (Polybutylenterephthalate) ●Resin ester polymer (carrier material) + pigment dispersion
C21918	PEG handle piece Luer pos.	●PBT (Polybutylenterephthalate) ●Resin ester polymer (carrier material) + pigment dispersion
C22125	LL neg. (cap)	●PBT (Polybutylenterephthalate) ●Resin ester polymer (carrier material) + pigment dispersion
C22126	Lug for PEG connection (cap)	●PA 6 (Polyamid) ●Resin ester polymer (carrier material) + pigment dispersion
C22127	Handle piece (cap)	●PBT (Polybutylenterephthalate) ●Resin ester polymer (carrier material) + pigment dispersion
J1040083		

Article code	Article name	Material
C21840	External fixation plate	●PDMS (Polydimethylesiloxane) (Silicone elastomer) containing vinyl groups, with Platinum catalyst and inhibitors ●Bariumsulphate ●PDMS (Polydimethylsiloxane) containing vinyl groups
C21841	Fixation clamp	●POM (Polyoxymethylene) or Acetal-Copolymer, base resin, additives ●Resin ester polymer (carrier material) + pigment dispersion
J1040084		
C21838	PEG handle piece Luer pos.	●PBT (Polybutylenterephthalate) ●Resin ester polymer (carrier material) + pigment dispersion
C21842	Connecting funnel	●TPU (thermoplastic Polyurethane), Polyether-based ●Carrier: resin-ester polymer with pigment dispersion
C22502	N-Heptan for Spectroscopy	●chemical
C22559	Glue	●ECA (Ethyl-2-Cyanoacrylate); Hydrochinon
J1040089		
C22200	Gastric retention disc	●PDMS (Polydimethylesiloxane) (Silicone elastomer) containing vinyl groups, with Platinum catalyst and inhibitors ●Bariumsulphate ●PDMS (Polydimethylsiloxane) containing Vinyl-groups
C22503	Ethanol absolute	●Ethanol absolute
C22559	Glue	●ECA (Ethyl-2-Cyanoacrylate); Hydrochinon
J1040066	tube injection CH 20	●TPU (thermoplastic Polyurethane), aromatic-based ●tube PUR (Polyurethane) 5,0x6,6mm ●tube printing ink: Ketone based colourant (cyclohexanone) ●tube printing ink t
J1040105		
C21843	Proximal cone	●TPU (thermoplastic Polyurethane), Polyester-based ●Resin ester polymer (carrier material) + pigment additive
C22559	Glue	●ECA (Ethyl-2-Cyanoacrylate); Hydrochinon
C22691	Fixation loop, OD 0.4 ± 0.01mm	●Core: PA 6.6 (Polyamide 6.6) ●Coating: PA 6 (Polyamide 6)
Packaging		
	Label Sticker 50 x 22 mm	●Paper and Adhesive

Article code	Article name	Material
	Adhesive label white 150 x 80 mm	• Upper basis paper: white Glassine paper B700-473, 55µm, 62 g/m² • Writing paper: satined article 1080, 74 g/m², 71µm The writing paper is coated from both sides with the adhesive • Adhesive: Permanent acrylic dispersion, 958 • Bottom basis paper: Glass
	Label sticker 140 x 175 mm	• Paper and Adhesive (acrylate dispersion)
	Natron Paper 1000 x 750 mm	• Color brown
	Imprinted envelope	• Printed paper (Velin 80 white, printed colour P 294)
	Instructions for use	• Paper
	Stretch film	• LLDPE (linear low density polyethylene)
	Stretch film foil plate	• LDPE (low density polyethylene)
	Cardboard 290 x 187 x 205 mm	• Corrugated cardboard
	Compound foil 30/70 µm	• PA (Polyamide)/PE (Polyethylene)
	Automat reel without print	• Paper
Sterilization		
	Indicator tape for EO sterilization	• Process indicator tape for Ethylene Oxide (EO) sterilization
	Sterilization gas	• Ethylene oxide