

Article code	7751562
Article name acc. DoC	Freka® Cone Adaptor, ENFit
Class acc. Directive 93/42/EEC	IIa

Article code	Article name	Material
C11503	Solvent	•THF (Tetrahydrofuran), also 1,4-Epoxy butane
C11515	Solvent	•MEK (Methylethylketone)
C11553	Solvent	•CHN (Cyclohexanone)
C22105	Urology connector	•PVC (Polyvinylchloride), colourant white, plasticizer DINP
C30823	ENFit female cap	•HDPE (high density Polyethylene) •colourant violet
C30824	Strap for ENFit cap, ring-Ø 10.5 mm	•LDPE (low density polyethylene), transparent
C30846	ENFit male, 21FR	•PVC (Polyvinylchloride) •colourant violet contains DOTE (Disubstituted Organic Tin); does not contain phthalate
C30959	Tube 5.1 x 7.1 mm	•PVC (Polyvinylchloride) with plasticizer DEHP, Colourant blue
C31525	Printing ink, blue	•Mixture (colourant and thinners)
C31526	Thinner for printing ink	•1-Methoxy-2-Propanol (2-Methoxymethylethoxy) propanol
Packaging		
	Roll Label 105 x 148 mm	•Paper (blank white) and Adhesive
	Natron Paper 1000 x 750 mm	•Color brown
	Stretch film	•LLDPE (linear low density polyethylene)
	Stretch film foil plate	•LDPE (low density polyethylene)
	Cardboard 249 x 148 x 150 mm	•Corrugated cardboard, color white
	Compound foil 30/70 µm	•PA (Polyamide)/PE (Polyethylene)
	Automat reel without print	•Paper
Sterilization		
	Indicator tape for EO sterilization	•Process indicator tape for Ethylene Oxide (EO) sterilization
	Sterilization gas	•100% Ethylene oxide