

Article code	7751937
Article name acc. DoC	GraviSet Varioline, ENFit
Class acc. MDR 2017/745	Is

Article code	Article name	Material
C10467	Roller clamp	●housing: PP (Polypropylene) plus colourant blue ●wheel: HDPE (high density polyethylene)
C11503	Solvent	●THF (Tetrahydrofuran), also 1,4-Epoxy butane
C31592	Tinopal Optical Brightener	●2,5-thiophenediylbis(5-tert-butyl-1,3-benzoxazole), color: Slight greenish yellow
C60977	PVC tube Ø3,5 x Ø5,0	●C90774: PVC (Polyvinylchloride) compound with plasticizer DEHT
J1010767		
C10330	Air ventilation plug	●C90330 ABS (Acrylonitrile Butadiene Styrene) ●C90405 Colourant white, ABS (Acrylonitrile Butadiene Styrene)-carrier
C10824	Bacteria filter	●glass microfiber filter
C10825	Bacteria filter band	●hydrophobic glass microfibre
C11120	Lid valve	●Silicone
C30032	Insert ENPlus	●C90351 HDPE (high density polyethylene) ●C90489 Colourant lilac
C30058	Body Varioline Adapter ENPlus	●C90383 LDPE (low density polyethylene) ●C90351 PE (Polyethylene) ●C90489 Colourant lilac ●C90286 lubricant demoulding aid
J1040192		
C30293	ENFit stopcock handle	●HDPE (high density polyethylene) with colourant
C30294	ENFit stopcock housing	●Copolyester (BPA-free)
C30295	Silicone oil SILFAR 350	●PDMS (Polydimethylsiloxane)
C30297	ENFit stopcock handle FK	●C90583 HDPE (high density polyethylene) ●C90584 colourbatch purple
C30298	ENFit stopcock housing FK	●C90598 Copolyester (BPA-free)
J1040205		
C30293	ENFit stopcock handle	●HDPE (high density polyethylene) with colourant
C30294	ENFit stopcock housing	●Copolyester (BPA-free)
C30295	Silicone oil SILFAR 350	●PDMS (Polydimethylsiloxane)
C30297	ENFit stopcock handle FK	●C90583 HDPE (high density polyethylene) ●C90584 colourbatch purple

Article code	Article name	Material
C30298	ENFit stopcock housing FK	●C90598 Copolyester (BPA-free)
J1040286		
C10333	Air ventilation plug	●C90735 ABS (Acrylonitrile-Butadiene-Styrene) ●C90210 Colourant white, ABS (Acrylonitrile-Butadiene-Styrene)-carrier
C10824	Bacteria filter	●glass microfiber filter
C11120	Lid valve	●Silicone
C30032	Insert ENPlus	●C90351 HDPE (high density polyethylene) ●C90489 Colourant lilac
C30058	Body Varioline Adapter ENPlus	●C90383 LDPE (low density polyethylene) ●C90351 PE (Polyethylene) ●C90489 Colourant lilac ●C90286 lubricant demoulding aid
J1040340		
C10897	Sealing for Cross spike	●LDPE (low density polyethylene) + EPE (foamed PE)
C11553	Solvent	●CHN (Cyclohexanone)
C30478	Union nut EN	●C90612 PP (Polypropylene) (50 %) ●C90452 Resin ester polymer (carrier material) + pigment dispersion
C30479	EN Drip chamber upper part	●C90735 ABS (Acrylonitrile-Butadiene-Styrene) ●C90734 Pigment ABS purple
C31180	Lower part of EN drip chamber for Ø3,0 x Ø5,0	●C90772: PVC (Polyvinylchloride) compound with plasticizer DINCH
J1040342		
C11499	Petroleum ether	●N/A
C11503	Solvent	●THF (Tetrahydrofuran), also 1,4-Epoxy butane
C11511	Lubricant	●Aminofunctional Siloxane
C30776	ENFit female 5.0 mm	●C90598: Copolyester (BPA-free) ●C90599: RAL 4005 (violet), Colour concentrate
C60977	PVC tube Ø3,5 x Ø5,0	●C90774: PVC (Polyvinylchloride) compound with plasticizer DEHT
Packaging		
C11020	Tubular foil for packaging	●Film 500m x 350 mm x 40 µm ●Mixture of different Polyolefines PE (Polyethylene), LDPE (Polyethylene, low density), transparent
C12184	Adhesive tape	●BOPP (biaxially oriented polypropylene) film and rubber adhesive

Article code	Article name	Material
C12704	Roll ENFit label 105×148.8 mm	●Paper (blank white with violette bar) and Adhesive
C12972	White paper for IFU printing	●Paper
C31138A	Carton 397 x 230 x 205 mm	●Paper
C31301	Compound foil 90µm	●Multilayer flat film PA (Polyamide) /adhesive layer / PE (Polyethylene)
C31579	PE-Band	●HDPE (high density polyethylene)
C64637	Peel paper GraviSets, ENFit ISO	●Paper
Sterilization		
C11520CN	Sterilization gas, Nanchang	●Ethylene oxide (90%) + Carbon dioxide (10%)
C11525	Indicator tape for EO sterilization	●Process indicator tape for Ethylene Oxide (EO) sterilization