

Article code	7751952
Article name acc. DoC	GraviSet EasyBag, universal, ENFit
Class acc. MDR 2017/745	Is

Article code	Article name	Material
C10467	Roller clamp	●housing: PP (Polypropylene) plus colourant blue ●wheel: HDPE (high density polyethylene)
C11503	Solvent	●THF (Tetrahydrofuran), also 1,4-Epoxy butane
C30832	ENFit male-ENLock Universal connector purple	●ABS (Acrylonitrile-Butadiene-Styrene) ●Colourant lilac
C30833	ENFit Male - LL Female Connector	●POM (Polyoxymethylene)
C30833A	ENFit male - LL female Connector	●PP (Polypropylene)
C30834	ENFit female oral / LL	●POM (Polyoxymethylene)
C60977	PVC tube Ø3,5 x Ø5,0	●C90774: PVC (Polyvinylchloride) compound with plasticizer DEHT
J1040192		
C30293	ENFit stopcock handle	●HDPE (high density polyethylene) with colourant
C30294	ENFit stopcock housing	●Copolyester (BPA-free)
C30295	Silicone oil SILFAR 350	●PDMS (Polydimethylsiloxane)
C30297	ENFit stopcock handle FK	●C90583 HDPE (high density polyethylene) ●C90584 colourbatch purple
C30298	ENFit stopcock housing FK	●C90598 Copolyester (BPA-free)
J1040205		
C30293	ENFit stopcock handle	●HDPE (high density polyethylene) with colourant
C30294	ENFit stopcock housing	●Copolyester (BPA-free)
C30295	Silicone oil SILFAR 350	●PDMS (Polydimethylsiloxane)
C30297	ENFit stopcock handle FK	●C90583 HDPE (high density polyethylene) ●C90584 colourbatch purple
C30298	ENFit stopcock housing FK	●C90598 Copolyester (BPA-free)
J1040220S		
C11499	Petroleum ether	●N/A
C11503	Solvent	●THF (Tetrahydrofuran), also 1,4-Epoxy butane
C11511	Lubricant	●Aminofunctional Siloxane
C30776	ENFit female 5.0 mm	●C90598: Copolyester (BPA-free) ●C90599: RAL 4005 (violet), Colour concentrate

Article code	Article name	Material
C30953	Tube 3.5 x 5.0 mm	●C90410 PVC (Polyvinylchloride) with plasticizer TOTM
J1040341		
C10594	Protective cap for ENPlus	●C90290 PP (Polypropylene)
C10897	Sealing for Cross spike	●LDPE (low density polyethylene) + EPE (foamed PE)
C11553	Solvent	●CHN (Cyclohexanone)
C30478	Union nut EN	●C90612 PP (Polypropylene) (50 %) ●C90452 Resin ester polymer (carrier material) + pigment dispersion
C30479	EN Drip chamber upper part	●C90735 ABS (Acrylonitrile-Butadiene-Styrene) ●C90734 Pigment ABS purple
C31180	Lower part of EN drip chamber for Ø3,0 x Ø5,0	●C90772: PVC (Polyvinylchloride) compound with plasticizer DINCH
J1040343		
C11499	Petroleum ether	●N/A
C11503	Solvent	●THF (Tetrahydrofuran), also 1,4-Epoxy butane
C11511	Lubricant	●Aminofunctional Siloxane
C30776	ENFit female 5.0 mm	●C90598: Copolyester (BPA-free) ●C90599: RAL 4005 (violet), Colour concentrate
C30832	ENFit male-ENLock Universal connector purple	●ABS (Acrylonitrile-Butadiene-Styrene) ●Colourant lilac
C60977	PVC tube Ø3,5 x Ø5,0	●C90774: PVC (Polyvinylchloride) compound with plasticizer DEHT
Packaging		
C11521	PE-Band 90 x 30 x 0.05 mm	●HDPE (high density polyethylene)
C12184	Adhesive tape	●BOPP (biaxially oriented polypropylene) film and rubber adhesive
C12704	Roll ENFit label 105×148.8 mm	●Paper (blank white with violette bar) and Adhesive
C12972	White paper for IFU printing	●Paper
C31301	Compound foil 90µm	●Multilayer flat film PA (Polyamide) /adhesive layer / PE (Polyethylene)
C51035A	Carton box	●Cardboard
C64634	Peel paper GraviSets, ENFit 100x225 mm	●Paper
Sterilization		
C11520CN	Sterilization gas, Nanchang	●Ethylene oxide (90%) + Carbon dioxide (10%)

Article code	Article name	Material
C11525	Indicator tape for EO sterilization	●Process indicator tape for Ethylene Oxide (EO) sterilization