

Article code	7751955
Article name acc. DoC	GraviSet EasyBag with cover, universal, ENFit
Class acc. MDR 2017/745	Is

Article code	Article name	Material
C10467	Roller clamp	●housing: PP (Polypropylene) plus colourant blue ●wheel: HDPE (high density polyethylene)
C10549	Closure cap LL male	●PP (Polypropylene) ●Masterbatch (white)
C10578	Protective cap	●PP (Polypropylene)
C11503	Solvent	●THF (Tetrahydrofuran), also 1,4-Epoxy butane,
C30833	ENFit Male - LL Female Connector	●POM (Polyoxymethylene)
C30834	ENFit female oral / LL	●POM (Polyoxymethylene)
C60977	PVC tube Ø3,5 x Ø5,0	●PVC (Polyvinylchloride) compound with plasticizer DEHT
J1040193		
	ENFit stopcock handle	●HDPE (high density polyethylene) with colourant
C30294	ENFit stopcock housing	●Copolyester (BPA-free)
C30295	Silicone oil SILFAR 350	●PDMS (Polydimethylsiloxane)
C30297	ENFit stopcock handle FK	●HDPE (high density polyethylene) ●colourbatch purple
C30298	ENFit stopcock housing FK	●Copolyester (BPA-free)
J1040206		
C30293	ENFit stopcock handle	●HDPE (high density polyethylene) with colourant
C30294	ENFit stopcock housing	●Copolyester (BPA-free)
C30295	Silicone oil SILFAR 350	●PDMS (Polydimethylsiloxane)
C30297	ENFit stopcock handle FK	●HDPE (high density polyethylene) ●colourbatch purple
C30298	ENFit stopcock housing FK	●Copolyester (BPA-free)
J1040341		
C10594	Protective cap for ENPlus	●PP (Polypropylene)
C10897	Sealing for Cross spike	●LDPE (low density polyethylene) + EPE (foamed PE)
C11553	Solvent	●CHN (Cyclohexanone)
C30478	Union nut EN	●PP (Polypropylene) (50 %) ●Resin ester polymer (carrier material) + pigment dispersion

Article code	Article name	Material
C30479	EN Drip chamber upper part	● ABS (Acrylonitrile-Butadiene-Styrene) ● C90734 Pigment ABS purple
	Lower part of EN drip chamber for Ø3,0 x Ø5,0	● PVC (Polyvinylchloride) compound with plasticizer DINCH
J1040343		
C11499	Petroleum ether	● N/A
C11503	Solvent	● THF (Tetrahydrofuran), also 1,4-Epoxy butane
	Lubricant	● Aminofunctional Siloxane
C30776	ENFit female 5.0 mm	● Copolyester (BPA-free) ● RAL 4005 (violet), Colour concentrate
C30832	ENFit male-ENLock Universal connector purple	● ABS (Acrylonitrile-Butadiene-Styrene) ● Colourant lilac
C60977	PVC tube Ø3,5 x Ø5,0	● PVC (Polyvinylchloride) compound with plasticizer DEHT
Packaging		
	PE-Band 90 x 30 x 0.05 mm	● HDPE (high density polyethylene)
	Adhesive tape	● BOPP (biaxially oriented polypropylene) film and rubber adhesive
	Roll ENFit label 105×148.8 mm	● Paper (blank white with violette bar) and Adhesive
	White paper for IFU printing	● Paper
	Compound foil 90µm	● Multilayer flat film PA (Polyamide) /adhesive layer / PE (Polyethylene)
	Carton box	● Cardboard
	top foil Gravi Sets ENFit 100x225mm	● Paper 60 g/m²
Sterilization		
	Sterilization gas, Nanchang	● Ethylene oxide (90%) + Carbon dioxide (10%)
	Indicator tape for EO sterilization	● Process indicator tape for Ethylene Oxide (EO) sterilization