

Article code	7752002
Article name acc. DoC	ProNeo Multifunction Syringe Cap
Class acc. Directive 93/42/EEC	Is

Article code	Article name	Material
C30289	Multifunction Syringe Cap	●syringe cap: HDPE (high density polyethylene), colourant - violet ●disc filter: ABS (Acrylonitrile Butadiene Styrene) + Nylon, colourant white
Packaging		
C11020	Tubular foil for packaging	●Film 500m x 350 mm x 40 µm ●Mixture of different Polyolefines PE (Polyethylene), LDPE (Polyethylene, low density), transparent
C11301	Compound foil 70 µm	●PA (Polyamide) / adhesive layer / PE (Polyethylene)
C12709	IFU ProNeo, A4	●copy paper 80g
C12768	Label sticker	●copy paper 80g
C31045	Stretch film	●LLDPE (linear low density polyethylene)
C31046	Stretch film foil plate	●LDPE (low density polyethylene)
	Carton 325 x 275 x 90 mm	●Cardboard
C64442	paper with single page and single colour printing; ProNeo Multifunction Syringe Cap	●Peel paper ●colour: black
Sterilization		
C11525	Indicator tape for EO sterilization	●Process indicator tape for Ethylene Oxide (EO) sterilization
C31500	Sterilization gas, Poland	●100% Ethylene oxide, liquefied grade one (GAZ 4-134) ●alternative: > 99.9% Ethylene oxide