

<b>Article code</b>	<b>7752003</b>
<b>Article name acc. DoC</b>	ProNeo MilkBottle Adapter
<b>Class acc. Directive 93/42/EEC</b>	Is

Article code	Article name	Material
C10894	Sealing ring	●PE (Polyethylene) foam
C30288	ProNeo Milk Bottle Adapter	●HDPE (high density polyethylene), colourant: C90353 blue lilac
C30759	Check valve with bacterial filter	<ul style="list-style-type: none"> <li>●Solvent: THF (Tetrahydrofuran)</li> <li>●Housing bottom part: PS (Polystyrene) + 2% Masterbatch PS (Polystyrene) lilac PS417925</li> <li>●Top part: PS (Polystyrene) + 2% Masterbatch PS (Polystyrene) universal white, UN MB0005</li> <li>●Filter holder: HDPE (high density Polyet</li> </ul>
C30823	ENFit female cap	<ul style="list-style-type: none"> <li>●HDPE (high density Polyethylene)</li> <li>●colourant violet</li> </ul>
C30897N	Strap for ENFit male - PEG	●N/A
<b>Packaging</b>		
C11081	Folding box 400 x 220 x 95 mm	●Corrugated cardboard
C12709	IFU ProNeo, A4	●copy paper 80g
	Label sticker	●copy paper 80g
C31045	Stretch film	●LLDPE (linear low density polyethylene)
C31046	Stretch film foil plate	●LDPE (low density polyethylene)
C31301	Compound foil 90µm	●Multilayer flat film PA (Polyamide) /adhesive layer / PE (Polyethylene)
C64419	paper with single page and single colour printing, ProNeo Milk Bottle Adaptor	<ul style="list-style-type: none"> <li>●Peel paper</li> <li>●colour: black</li> </ul>
<b>Sterilization</b>		
C11525	Indicator tape for EO sterilization	●Process indicator tape for Ethylene Oxide (EO) sterilization
C31500	Sterilization gas, Poland	<ul style="list-style-type: none"> <li>●100% Ethylene oxide, liquefied grade one (GAZ 4-134)</li> <li>●alternative: &gt; 99.9% Ethylene oxide</li> </ul>