

Article code	7752003
Article name acc. DoC	ProNeo MilkBottle Adapter
Class acc. Directive 93/42/EEC	Is

Article code	Article name	Material
C10894	Sealing ring	●PE (Polyethylene) foam
C30288	ProNeo Milk Bottle Adapter	●HDPE (high density polyethylene), colourant: blue lilac
C30759	Check valve with bacterial filter	●Solvent: THF (Tetrahydrofuran) ●Housing bottom part: PS (Polystyrene) + Masterbatch PS (Polystyrene) lilac PS417925 ●Top part: PS (Polystyrene) + Masterbatch PS (Polystyrene) universal white, UN MB0005 ●Filter holder: HDPE (high density Polyet
C30823	ENFit female cap	●HDPE (high density Polyethylene) ●colourant violet
C30897N	Strap for ENFit male - PEG	●N/A
Packaging		
	Folding box 400 x 220 x 95 mm	●Corrugated cardboard
	IFU ProNeo, A4	●copy paper 80g
	Label sticker	●copy paper 80g
	Stretch film	●LLDPE (linear low density polyethylene)
	Stretch film foil plate	●LDPE (low density polyethylene)
	Compound foil 90µm	●Multilayer flat film PA (Polyamide) /adhesive layer / PE (Polyethylene)
	paper with single page and single colour printing, ProNeo Milk Bottle Adaptor	●Peel paper ●colour: black
Sterilization		
	Sterilization gas	●Ethylene oxide
	Indicator tape for EO sterilization	●Process indicator tape for Ethylene Oxide (EO) sterilization