

Article code	7752004
Article name acc. DoC	ProNeo Medication Adaptor
Class acc. Directive 93/42/EEC	Is

Article code	Article name	Material
C30823	ENFit female cap	● HDPE (high density Polyethylene) ● colourant violet
C30897N	Strap for ENFit male - PEG	● N/A
J1020086		
C10581	Protective Cap	● PP (Polypropylene)
C10884	Air filter 0.2 µm	● ePTFE (expanded Polytetrafluor-ethylene) ● carrier: Polyester fleece
C21812	Coverage for air filter	● ABS (Acrylonitrile Butadiene Styrene) ● Colorbatch RAL4005 Terlux 2802HD blue lilac MABS (Methymethacrylate Acrylonitrile Butadiene Styrene)-basis + pigment
C21813	Filter coverage for Medication Adapter	● ABS (Acrylonitrile Butadiene Styrene) ● MABS (Methymethacrylate Acrylonitrile Butadiene Styrene) RAL 4005
C30286	ProNeo Medication Adapter Housing	● ABS (Acrylonitrile Butadiene Styrene) ● Colorbatch RAL4005 Terlux 2802HD blue lilac MABS-basis + pigment
Packaging		
	Tubular foil for packaging	● Film 500m x 350 mm x 40 µm ● Mixture of different Polyolefines PE (Polyethylene), LDPE (Polyethylene, low density), transparent
	Folding box 400 x 220 x 95 mm	● Cardboard, color white outside and brown inside
	Compound foil 70 µm	● PA (Polyamide) / adhesive layer / PE (Polyethylene)
	IFU ProNeo, A4	● copy paper 80g
	Label sticker	● copy paper 80g
	Stretch film	● LLDPE (linear low density polyethylene)
	Stretch film foil plate	● LDPE (low density polyethylene)
	paper with single page and single colour printing; ProNeo Medication Adaptor	● Peel paper ● colour: black
Sterilization		
	Indicator tape for EO sterilization	● Process indicator tape for Ethylene Oxide (EO) sterilization
	Sterilization gas	● 100% Ethylene oxide