

<b>Article code</b>	7752013
<b>Article name acc. DoC</b>	ENFit Sampling Straw, 150 mm
<b>Class acc. Directive 93/42/EEC</b>	Is

<b>Article name</b>	<b>Material</b>
Straw	PUR
Hub Connector	MABS
Sheath	Polypropylene
Packaging	
Primary Packaging	Porous Material – Paper & PET/Pe Composite
Secondary Packaging	Hard Paper
Sterilization	
Sterilization gas	●100% Ethylene oxide