

Article code	7752015
Article name acc. DoC	ENFit VarioBottle Adapter
Class acc. Directive 93/42/EEC	Is

Article name	Material
Cap	HDPE + colorant
Strap	LDPE
Connector	ABS + colorant
Multisteppped bottle adapter	MARFRAN MED 65A ST + colorant
Bush	Styrene-Butadiene Copolymer + colorant
Packaging	
Primary Packaging – dust proof paper and sealed blister	Paper, PE-PA
Secondary Packaging	Hard Paper
Sterilization	
Sterilization gas	100% Ethylene oxide