

Article code	7752030
Article name acc. DoC	ProNeo Feeding Tube FR 8, 50 cm
Class acc. Directive 93/42/EEC	IIa

Article code	Article name	Material
C11503	Solvent	●THF (Tetrahydrofuran), also 1,4-Epoxy butane
C21948	Tube tip	●PUR (Polyurethane) with 10% BaSO ₄
C22501	Aceton - for analysis	●Acetone (2-Propanone, Dimethyl ketone)
C22572	Lubricant	●PDMS (Polydimethylsiloxane)
C22720	tube plaster, heart 75 mm x 33 mm x 0.32 mm	●one layer of polyester fleece, non-woven polyacrylateadhesive siliconized paper liner, white
C22973	Guide wire ENFit ProNeo 505mm	●C22639 Guide wire: Stainless steel 1.4310 (DIN EN 10270-3) wire ●C30885 ENFit wire grip: MABS (Methylmethacrylate Acrylnitrile Butadiene Styrene) + Colorbatch ●C30819 Adhesive
C30280	ProNeo Connector Male FR 8	●MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene) ●colourant: blue lilac MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene-basis + pigment)
C30352	ProNeo Soft Cap, white	●POE (Polyolefine elastomers) ●lubricant copolymer of ethylene + octene ●colourant - white
C30981	Tube 1.9 x 2.85 mm	●PUR (Polyurethane) and PUR (Polyurethane) with 60% BaSO ₄ + Colourant (white)
C31534	Tube printing ink	●Butanone ●Solvent Black 29 ●Ethanol ●Cellulose Nitrate ●Isopropyle Acetate ●2-Naphtole
C31535	Tube printing ink thinner	●Butanone ●Acetone
Packaging		
	Edge protection	●corrugated cardboard
	Compound foil 70 µm	●PA (Polyamide) / adhesive layer / PE (Polyethylene)
	Peel paper	●Paper 60 g/m ²
	Label sticker	●copy paper 80g
	IFU ProNeo, A4	●copy paper 80g
	Natron Paper 1000 x 750 mm	●Color brown
	Stretch film	●LLDPE (linear low density polyethylene)

Article code	Article name	Material
	Stretch film foil plate	●LDPE (low density polyethylene)
	Cardboard 249 x 148 x 100 mm	●Corrugated cardboard, color white
	Black ink	●Printing ink
Sterilization		
	Indicator tape for EO sterilization	●Process indicator tape for Ethylene Oxide (EO) sterilization
	Sterilization gas, Poland	●Ethylene oxide