

|                                       |                           |
|---------------------------------------|---------------------------|
| <b>Article code</b>                   | 7755643                   |
| <b>Article name acc. DoC</b>          | Freka® PEG CH/FR15, ENFit |
| <b>Class acc. Directive 93/42/EEC</b> | IIb                       |

| Article code | Article name                     | Material                                                                                                                             |
|--------------|----------------------------------|--------------------------------------------------------------------------------------------------------------------------------------|
| C11515       | Solvent                          | ●MEK (Methylethylketone)                                                                                                             |
| C11553       | Solvent                          | ●CHN (Cyclohexanone)                                                                                                                 |
| C21832       | Spool and thread                 | ● MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene), Colourant (blue)<br>● PA 6 / PA 6.6                                     |
| C21845       | Thread introducer device         | ● PP (Polypropylene)<br>● Resin ester polymer (carrier material) + pigment dispersion                                                |
| C21941       | Tube clamp with centering, white | ●POM (Polyoxymethylene)<br>●Color Carrier: resin-ester polymer with pigment dispersion                                               |
| C22606       | Tray cover sheet                 | ●coloured PETG (Polyethyleneterephthalate Glycol), blue tinted, mould release: silicone additive                                     |
| C22609       | PEG tray 200 x 155 mm            | ●coloured PETG (Polyethyleneterephthalate Glycol), blue tinted, mould release: silicone additive                                     |
| C31526       | Thinner for printing ink         | ●1-Methoxy-2-Propanol (2-Metholymethylethoxy) propanol                                                                               |
| C31527       | Printing ink, black              | ●Stamping ink                                                                                                                        |
| C31586       | PenBlade-11P Safety Scalpel      | ●blade: stainless steel<br>●body: ABS (Acrylonitrile Butadiene Styrene) resin, colourant: orange + white<br>●spring: SUS 304 Polymer |
| C22059       | Membrane                         | ● LSR (Liquid Silicone Rubber), biomedical grade                                                                                     |
| C22061       | Valve upper part                 | ● PC (Polycarbonate)<br>● Dispersion of pigments in a universal carrier (resin ester-polymer)                                        |
| C22062       | Valve lower part                 | ● PC (Polycarbonate)<br>● Masterbatch: Colourant, white, 3%                                                                          |
| C22598       | Spring                           | ●Stainless steel 1.4401                                                                                                              |
| C22638       | Ethanol 70%                      | ●Ethyl alcohol 70%                                                                                                                   |
| J1040029     |                                  |                                                                                                                                      |
| C11515       | Solvent                          | ●MEK (Methylethylketone)                                                                                                             |
| C11553       | Solvent                          | ●CHN (Cyclohexanone)                                                                                                                 |
| C21803       | Blood stop plug                  | ●MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene)                                                                           |
| C22081       | Valve, lower part                | ●PUR (Polyurethane) with masterbatch green and release additive                                                                      |
| C22082       | Cannula valve upper part         | ●PUR (Polyurethane)<br>●Colourant                                                                                                    |
| C22084       | Valve opening piece              | ●PP (Polypropylene)                                                                                                                  |
| C22085       | Protective cap                   | ●PE (Polyethylene)                                                                                                                   |

| Article code | Article name                                       | Material                                                                                                                                                                                                   |
|--------------|----------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| C22558       | Glue                                               | ●ECA (Ethyl-2-Cyanoacrylate); Hydrochinon (0.01- < 0.1% w/w)                                                                                                                                               |
| C22573       | Lubricant                                          | ●Siloxane                                                                                                                                                                                                  |
| C22588       | Cannula                                            | ●needle holder: PC (polycarbonate based on BPA bisphenole A)<br>●needle: stainless steel<br>●glue Polyurethane                                                                                             |
| C22665       | Capillary 8/80                                     | ●Fluorinated Ethylene Propylene (FEP);<br>●Stripes: Pantone 542 blue;<br>●X-ray contrast agent: 24% BaSO <sub>4</sub> (pre-mixed in main material)                                                         |
| C22937S      | Membrane                                           | ●LSR (Liquid silicone rubber)                                                                                                                                                                              |
| C31757       | PTFE<br>(Polytetrafluoroethylene)<br>Etching Agent | ●Etching Agent for PTFE                                                                                                                                                                                    |
| J1040088     |                                                    |                                                                                                                                                                                                            |
| C22201       | Gastric retention disc                             | ● PDMS (Polydimethylsiloxane) (Silicone elastomer) containing vinyl groups, with Platinum catalyst and inhibitors<br>● Bariumsulphate,<br>● PDMS (Polydimethylsiloxane) containing Vinyl-groups            |
| C22503       | Ethanol absolute                                   | ●Ethanol absolute                                                                                                                                                                                          |
| C22559       | Glue                                               | ●ECA (Ethyl-2-Cyanoacrylate); Hydrochinon (0.01- < 0.1% w/w)                                                                                                                                               |
| J1040063     | tube injection CH 15                               | ● TPU (thermoplastic Polyurethane), aromatic-based<br>●tube PUR (Polyurethane) 3,65x4,8mm<br>● tube printing ink: Ketone based colourant (cyclohexanone (5-10%), butanone (70-80%))<br>● tube printing ink |
| J1040090     |                                                    |                                                                                                                                                                                                            |
| C21891       | External fixation plate                            | ● PDMS (Polydimethylsiloxane) (Silicone elastomer) containing vinyl groups, with Platinum catalyst and inhibitors<br>● Bariumsulphate<br>● PDMS (Polydimethylsiloxane) containing vinyl groups             |
| C21893       | Fixation clamp                                     | ● POM (Polyoxymethylene) or Acetal-Copolymer, base resin (≥ 99%);, additives (< 1%)<br>● Resin ester polymer (carrier material) + pigment dispersion                                                       |
| J1040104     |                                                    |                                                                                                                                                                                                            |
| C21872       | Proximal cone                                      | ● TPU (thermoplastic Polyurethane), Polyester-based<br>● Carrier material: resin-ester polymer with pigment additive                                                                                       |
| C22559       | Glue                                               | ●ECA (Ethyl-2-Cyanoacrylate); Hydrochinon (0.01- < 0.1% w/w)                                                                                                                                               |
| C22691       | Fixation loop, OD 0.4 ± 0.01mm                     | ●Core: PA 6.6 (Polyamide 6.6)<br>●Coating: PA 6 (Polyamide 6)                                                                                                                                              |

| Article code  | Article name                            | Material                                                                                                                                                                               |
|---------------|-----------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| J1040196      |                                         |                                                                                                                                                                                        |
| C30823        | ENFit female cap                        | <ul style="list-style-type: none"> <li>● HDPE (high density Polyethylene)</li> <li>● colourant violet</li> </ul>                                                                       |
| C30892        | ENFit male - PEG connector              | <ul style="list-style-type: none"> <li>● PBT (Polybutylenterephthalate)</li> <li>● Resin ester polymer (carrier material)</li> </ul>                                                   |
| C30897        | Strap for ENFit male connector          | <ul style="list-style-type: none"> <li>● LDPE (Low density polyethylene),</li> </ul>                                                                                                   |
| J1040207      |                                         |                                                                                                                                                                                        |
| C21857        | Fixation aid PEG                        | <ul style="list-style-type: none"> <li>● PBT (Polybutylenterephthalate)</li> <li>● Resin ester polymer (carrier material) + pigment dispersion</li> </ul>                              |
| C21859        | Fixation ring                           | <ul style="list-style-type: none"> <li>● TPE (Thermoplastic Elastomer)</li> <li>● color: white</li> </ul>                                                                              |
| C21861        | Fixation screw                          | <ul style="list-style-type: none"> <li>● PBT (Polybutyleneterephthalate)</li> <li>● Resin ester polymer (carrier material) + pigment dispersion</li> </ul>                             |
| Packaging     |                                         |                                                                                                                                                                                        |
|               | Adhesive tape paper                     | <ul style="list-style-type: none"> <li>● gummed paper tapes</li> <li>● Polyeste filament yarn for strengthening</li> </ul>                                                             |
|               | Label sticker, size: 97 x 67.7 mm       | <ul style="list-style-type: none"> <li>● Paper 70 g / qm</li> <li>● adhesive (Acrylate-dispersion)</li> </ul>                                                                          |
|               | Tamper-evident safety seal              | <ul style="list-style-type: none"> <li>● PP (Polypropylene) -Film Transparent- 60 µm Top Coat</li> <li>● Adhesive: Permanent -958</li> <li>● Base: Glassine paper B 700-476</li> </ul> |
|               | Label Sticker 50 x 22 mm                | <ul style="list-style-type: none"> <li>● Paper and Adhesive</li> </ul>                                                                                                                 |
|               | Label sticker 140 x 175 mm              | <ul style="list-style-type: none"> <li>● Paper and Adhesive (acrylate dispersion)</li> </ul>                                                                                           |
|               | Instructions for use                    | <ul style="list-style-type: none"> <li>● Paper</li> </ul>                                                                                                                              |
|               | Stretch film                            | <ul style="list-style-type: none"> <li>● LLDPE (linear low density polyethylene)</li> </ul>                                                                                            |
|               | Stretch film foil plate                 | <ul style="list-style-type: none"> <li>● LDPE (low density polyethylene)</li> </ul>                                                                                                    |
|               | Cardboard 324 x 279 x 188 mm            | <ul style="list-style-type: none"> <li>● Corrugated cardboard</li> </ul>                                                                                                               |
|               | Compound foil 30/70 µm                  | <ul style="list-style-type: none"> <li>● PA (Polyamide)/PE (Polyethylene)</li> </ul>                                                                                                   |
|               | Folded carton box, brown, 280x180x45 mm | <ul style="list-style-type: none"> <li>● Material: Testliner</li> <li>● Wave: E</li> </ul>                                                                                             |
|               | Automat reel without print              | <ul style="list-style-type: none"> <li>● Paper</li> </ul>                                                                                                                              |
| Sterilization |                                         |                                                                                                                                                                                        |
|               | Indicator tape for EO sterilization     | <ul style="list-style-type: none"> <li>● Process indicator tape for Ethylene Oxide (EO) sterilization</li> </ul>                                                                       |
|               | Sterilization gas                       | <ul style="list-style-type: none"> <li>● 100% Ethylene oxide</li> </ul>                                                                                                                |