

Article code	7755647 / M90800310
Article name acc. DoC	Freka Intestinal Tube FR 12 for Freka® PEG FR 20, ENFit
Class acc. Directive 93/42/EEC	IIB

Article code	Article name	Material
C11515	Solvent	●MEK (Methylethylketone)
C11553	Solvent	●CHN (Cyclohexanone)
C21839	Fixation ring	●SEBS Compound (Styrene-Ethylene / Butylene-Styrene)
C21857	Fixation aid PEG	●PBT (Polybutylenterephthalate) ●Resin ester polymer (carrier material) + pigment dispersion
C21867	Closure cap Luer Lock pos.	●TPU (thermoplastic Polyurethane), Polyether-based ●Preparation of colour agent in a universal carrier
C21937	End piece (flexible tip)	●PUR (Polyurethane) with 10% BaSO ₄ ●Lubricant (Amidwachs) ●Slip additive Colorsip 103 100% carboxylic acid ester
C22503	Ethanol absolute	●Ethanol absolute
C22537	Guide wire, 1.14 x 1460 mm	●Stainless steel, coating: PTFE (Polytetrafluoroethylene)
C22572	Lubricant	●PDMS (Polydimethylsiloxane)
C30653	Luer Lock pos.	●TPU (thermoplastic Polyurethane), Polyester-based ●Resin ester polymer (carrier material) + pigment additive
C30812	Fixation screw	●PBT (Polybutylenterephthalate) ●Resin ester polymer (carrier material) + pigment dispersion
C30972	Tube 6.1 mm as part of glue mixture	●PUR (Polyurethane) and PUR + 60% BaSO ₄ (one x-ray contrast strip)
C30989	Tube 2.8 x 4.05 x 1200 mm	●PUR (Polyurethane) and PUR + 60% BaSO ₄ (one x-ray contrast strip)
C31526	Thinner for printing ink	●1-Methoxy-2-Propanol (2-Methoxymethylethoxy) propanol
C31527	Printing ink, black	●Stamping ink
C31534	Tube printing ink	●Butanone 70-80% C.I ●Solvent Black <10% ●Ethanol 5-10% ●Cellulose Nitrate 5-10% ●Isopropyl Acetate 5% ●2-Naphthol <0.2%
C31535	Tube printing ink thinner	●Butanone 96-98% ●Acetone 2-4%
J1040110		

Article code	Article name	Material
C21831	Interior fixation screw	●PBT (Polybutylenterephthalate) ●Resin ester polymer (carrier material) + pigment dispersion
C21849	Adaptor exterior part	●POM (Polyoxymethylene) or Acetal-Copolymer, base resin (≥ 99%), additives (< 1%) ●Carrier: resin-ester polymer with pigment dispersion
C21907	Fixation screw	●PBT (Polybutylenterephthalate) ●Resin ester polymer (carrier material) + pigment dispersion
C21936	Fixation ring	●TPE (Thermoplastic Elastomer) ● color: white
C21938	Adaptor interior part	●PBT (Polybutylenterephthalate) ●Carrier: resin-ester polymer with pigment dispersion
C22054	O-ring	●Siloxane MVQ-Silicone (Methyl-Vinyl-Silicone rubber)
C22559	Glue	●ECA (Ethyl-2-Cyanoacrylate)
J1040135		
C21811	Connecting funnel	●TPU (thermoplastic Polyurethane), Polyether-based ●Carrier: resin-ester polymer with pigment dispersion ●Lubricant (Amid wax)
C21833	Sealing body for click adaptor	●PBT (Polybutylenterephthalate) ●Resin ester polymer (carrier material) + pigment dispersion
C22559	Glue	●ECA (Ethyl-2-Cyanoacrylate)
C30808	ENLock Y-piece	●PBT (Polybutylenterephthalate) ●Resin ester polymer (carrier material) + pigment dispersion
J1040201		
C22562	Mini tube 3.2 x 0.2 x 11.5 mm	●Stainless steel
C22596	Glue	●Epoxy Resin
C22597	Glue hardener	●Modified Polyamine
C30823	ENFit female cap	●HDPE (high density Polyethylene) ●colourant violet
C30896	ENFit male - PEG connector	●PBT (Polybutylenterephthalate) ●Resin ester polymer (carrier material)
C30897	Strap for ENFit male connector	●PE (Polyethylene)
J1060209		
C21988	Adaptor with Luer Lock neg.	●PA 6 (Polyamide)
C21989	Seal ring 3.1 x 6.5 x 1.6 mm	●PE (Polyethylene)
C21990	Adaptor insertion	●SEBS (Styrene-ethylene/ butylene-styrene)
C21991	Y-piece	●PC (Polycarbonate)

Article code	Article name	Material
Packaging		
	Tamper-evident safety seal	●PP (Polypropylene) -Film Transparent- 60 µm Top Coat ●Adhesive: Permanent -958 ●Base: Glassine paper B 700-476
	Label Sticker 50 x 22 mm	●Paper and Adhesive
	Label sticker 140 x 175 mm	●Paper and Adhesive
	Instructions for use	●Paper
	Stretch film	●LLDPE (linear low density polyethylene)
	Stretch film foil plate	●LDPE (low density polyethylene)
	Cardboard 290 x 187 x 205 mm	●Corrugated cardboard, color brown
	Compound foil 30/70 µm	●PA (Polyamide)/PE (Polyethylene)
	Folded carton box, brown, 270x150x30 mm	● Material: Testliner ● Wave: E
	Automat reel without print	●Paper
Sterilization		
	Indicator tape for EO sterilization	●Process indicator tape for Ethylene Oxide (EO) sterilization
	Sterilization gas	●100% Ethylene oxide