

Article code	7755648 (M90800311)
Article name acc. DoC	Freka GastroTube FR 15, ENFit
Class acc. Directive 93/42/EEC	IIb

Article code	Article name	Material
C21887	Fixation plate	●PDMS (Polydimethylsiloxane) with Vinyl Groups
C21888	Fixation clamp	●POM (Polyoxymethylene) or Acetal-Copolymer, base resin (≥ 99%), additives (< 1%)
C21935	Tube clamp	●POM (Polyoxymethylene) or Acetal-Copolymer, base resin (≥ 99%), additives (< 1%) ●Resin ester polymer (carrier material) + pigment dispersion
C22014	Guide wire with handhold gastro tube	●MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene) ●Masterbatch ivory
C22572	Lubricant	●PDMS (Polydimethylsiloxane)
C22595	Ring Ø 5 x thickness 0.1 x height 1 mm	●Silver, medical grade
C22638	Ethanol 70%	●Ethyl alcohol 70%
C22668	Syringe 10 ml	●Piston: PE (Polyethylene), PP (Polypropylene) ●Cylinder: PP (Polypropylene) + white colourant
C22999	ENFit Replacement Set Body	●LSR (Liquid Silicone Rubber), biomedical grade with one stripe of platinum cure silicone with 11% radioopaque ingredients
C23157	Glue	●Silicone medical grade
C30419	ENFit male insert for GastroTube and SIK-tube funnel	●MABS (Methylmethacrylate-Acrylonitrile-Butadiene-Styrene), ● Masterbatch violett (RAL4005)
J1040017		
C22059	Membrane	●LSR (Liquid Silicone Rubber), biomedical grade
C22061	Valve upper part	●PC (Polycarbonate) ●Dispersion of pigments in a universal carrier (resin ester-polymer)
C22062	Valve lower part	●PC (Polycarbonate)
C22598	Spring	●Stainless steel 1.4401
C22638	Ethanol 70%	●Ethyl alcohol 70%
J1040027		
C22174	MCT-oil and vials	●Oil: medium-chain triglycerides - caproic acid 0-2%, caprylic acid 50-65%, capric acid 35-45%, lauric acid 0-3% ●Vials: LDPE (low density polyethylene)
C22636	Seal bags 145 x 60 mm	●PES (Polyester) ●ALU (Aluminium) foil ●PE (Polyethylene)

Article code	Article name	Material
Packaging		
	Label Sticker 50 x 22 mm	●Paper and Adhesive
	Adhesive label white 150 x 80 mm	● Upper basis paper: white Glassine paper B700-473, 55µm, 62 g/m² ● Writing paper: satined article 1080, 74 g/m², 71µm The writing paper is coated from both sides with the adhesive ● Adhesive: Permanent acrylic dispersion, 958 ● Bottom basis paper: Glass
	Label sticker 140 x 175 mm	●Paper and Adhesive
	Instructions for use for 7755648	●Paper
	Self-adhesive tape, white	● PP (Polypropylene) ● adhesive: natural rubber
	Natron Paper 1000 x 750 mm	●Color brown
	Stretch film	●LLDPE (linear low density polyethylene)
	Stretch film foil plate	●LDPE (low density polyethylene)
	Folding box 782 x 595 x 360 mm	●Corrugated cardboard
	Cardboard 253 x 150 x 85 mm	●Corrugated cardboard, color white
	Compound foil 30/70 µm	●PA (Polyamide)/PE (Polyethylene)
	Folded carton box, brown, 270x150x30 mm	● Material: Testliner ● Wave: E
	Automat reel without print	●Paper
Sterilization		
	Indicator tape for EO sterilization	●Process indicator tape for Ethylene Oxide (EO) sterilization
	Sterilization gas	●100% Ethylene oxide