

Article code	7755693
Article name acc. DoC	Giving Set Extension 100 cm, ENFit
Class acc. Directive 93/42/EEC	IIa

Article code	Article name	Material
C10463	Roberts clamp, white	●ABS (Acrylonitrile Butadiene Styrene)
C10942 C10942A	Tube 3.0 x 4.1 mm	●Polyether Polyurethane
C11499	Petroleum ether	●N/A
C11503	Solvent	●THF (Tetrahydrofuran), also 1,4-Epoxy butane
C11511	Lubricant	●Aminofunctional Siloxane
C30811	ENFit female 4.1 mm	●Tritan Copolyester, transparent ●RAL 4005 (violet), Colour concentrate for Tritan MX731
C30842	ENFit male, 12FR	●PET (Polyethyleneterephthalat) (Copolyester) ●colourant violet
Packaging		
	Compound foil 70 µm	●PA (Polyamide) / adhesive layer / PE (Polyethylene)
	Peel paper	●Paper 60 g/m2
	Natron Paper 1000 x 750 mm	●Color brown
	Label sticker 160 x 115 mm	●Paper
	Stretch film	●LLDPE (linear low density polyethylene)
	Stretch film foil plate	●LDPE (low density polyethylene)
	Cardboard 290 x 187 x 205 mm	●Corrugated cardboard
	Black ink	●Printing ink
	Blue ink	●Printing ink
Sterilization		
	Indicator tape for EO sterilization	●Process indicator tape for Ethylene Oxide (EO) sterilization
	Sterilization gas	●100% Ethylene oxide