

Article code	7755695
Article name acc. DoC	Freka® Universal Funnel Adaptor, ENFit
Class acc. MDR 2017/745	I

Article code	Article name	Material
C11503	Solvent	•THF (Tetrahydrofuran), also 1,4-Epoxy butane
C22023	Universal funnel adaptor	•TPU (thermoplastic Polyurethane), Polyether-based •Resin ester polymer (carrier material) + pigment dispersion •Lubricant (Amid wax)
C30835	ENFit Female - LL Male Connector	•PP (Polypropylene)
Packaging		
C12704	Roll ENFit label 105×148.8 mm	•Paper (blank white with violette bar) and Adhesive
C12778	Additional information	Paper
C22645	Natron Paper 1000 x 750 mm	•Color brown
C31045	Stretch film	•LLDPE (linear low density polyethylene)
	Stretch film foil plate	•LDPE (low density polyethylene)
C31064A	Cardboard 290 x 187 x 205 mm	•Corrugated cardboard
C31094	Zip lock bag	•PE (Polyethylene)