

Article code	7840152AP
Article name acc. DoC	Freka® Sonde CH 10, 120 cm, ENLock
Class acc. Directive 93/42/EEC	IIa

Article code	Article name	Material
C11503	Solvent	●THF (Tetrahydrofuran), also 1,4-Epoxy butane
C22501	Aceton - for analysis	●Acetone (2-Propanone, Dimethyl ketone)
C22572	Lubricant	●PDMS (Polydimethylsiloxane)
C22704	Fixation tape	●polyester fleece non-woven ●polyacrylate adhesive ●siliconized paper liner
C30686	Cap	●POE (Polyolefine elastomers) ●lubricant copolymer of ethylene + octene ●carrier material + pigment dispersion
C30728	Tube Tip FR 10	●PUR (Polyurethane) with 10% BaSO4 ●Lubricant (Amidwachs) ●Slip additive Colorslip 103 100 % carboxylic acid ester
C30875	ENLock Universal Y Funnel FR 10	● PUR (Polyurethane) ● glue: CHN (Cyclohexanone)
C30906	Tube 2.4 x 3.4 mm	●PUR (Polyurethane) + PUR with 60% BaSO4 + Colourant (white)
C31527	Printing ink, black	●Stamping ink
C31534	Tube printing ink	●Butanone ●Solvent Black ●Ethanol ●Cellulose Nitrate ●Isopropyle Acetate ●2-Naphtole
C31535	Tube printing ink thinner	●Butanone ●Acetone
J1040120		
C22633	Mandrin 1.37 x 1205 mm	●Stainless steel
C30814	Handle Piece	●MABS (Methylmethacrylate Acrylnitrile Butadiene Styrene) ●Resin ester polymer (carrier material) + pigment dispersion
C30819	Glue	●superglue based on Alkoxyethyl cyanoacrylate
Packaging		
	White paper for IFU printing	●100% wood pulp
	Label sticker 145 x 94 mm	●80 g/sqm adhesive paper (adhesive based on acrylic)
	Natron Paper 1000 x 750 mm	●Color brown
	Stretch film	●LLDPE (linear low density polyethylene)
	Stretch film foil plate	●LDPE (low density polyethylene)

Article code	Article name	Material
	Cardboard 249 x 148 x 100 mm	●Corrugated cardboard, color brown ●Printing FK Logo: Blue- Pantone 300U
	Compound foil 30/70 µm	●PA (Polyamide)/PE (Polyethylene)
	Thinner for printing ink	●1-Methoxy-2-Propanol (2-Metholymethylethoxy) propanol
	Paper	●Paper
Sterilization		
	Indicator tape for EO sterilization	●Process indicator tape for Ethylene Oxide (EO) sterilization
	Sterilization gas	●Ethylene oxide