

Article code	7901111AP
Article name acc. DoC	Freka PEG Set gastric FR 15
Class acc. Directive 93/42/EEC	I Ib

Article code	Article name	Material
C21832	Spool and thread	● MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene), Colourant (blue) ● PA 6 / PA 6.6 (Polyamid)
C21845	Thread introducer device	● PP (Polypropylene) ● Resin ester polymer (carrier material) + pigment dispersion
C21881	Connecting funnel with Luer	● TPU (thermoplastic Polyurethane), Polyether-based ● Resin ester polymer (carrier material) + pigment dispersion ● Lubricant (Amid wax)
C21941	Tube clamp with centering, white	● POM (Polyoxymethylene) ● Color Carrier: resin-ester polymer with pigment dispersion
C22606	Tray cover sheet	● coloured PETG (Polyethyleneterephthalate Glycol), blue tinted, mould release: silicone additive
C22609	PEG tray 200 x 155 mm	● coloured PETG (Polyethyleneterephthalate Glycol), blue tinted, mould release: silicone additive
C22663	Scalpel	● Blade: stainless steel ● Handle: HIPS (high impact polystyrene) ● Guard: LDPE (low density polyethylene)
C22918	locking screw FR 15	● PBT (Polybutyleneterephthalate), color batch: universal Masterbatch ● PBT (Polybutyleneterephthalate) ● TPE (thermoplastic elastomere)
C30976	Tube 3.65 x 4.8 x 350 ± 5 mm	● PUR (Polyurethane) and PUR + 60% BaSO₄ (one x-ray contrast strip)
C31531	Tube printing ink	● Ketone based colourant (cyclohexanone , butanone)
C31532	Tube printing ink thinner	● Ketone based additive (Butanone = Methyllethylketone)
J1040029		
C11515	Solvent	● MEK (Methyllethylketone)
C11553	Solvent	● CHN (Cyclohexanone)
C21803	Blood stop plug	● MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene)
C22081	Valve, lower part	● PUR (Polyurethane) with masterbatch green and release additive
C22082	Cannula valve upper part	● PUR (Polyurethane) ● Colourant

Article code	Article name	Material
C22084	Valve opening piece	●PP (Polypropylene)
C22085	Protective cap	●PE (Polyethylene)
C22558	Glue	●ECA (Ethyl-2-Cyanoacrylate); Hydrochinon
C22573	Lubricant	●Siloxane
C22588	Cannula	●needle holder: PC (polycarbonate based on BPA bisphenole A) ●needle: stainless steel ●glue Polyurethane
C22665	Capillary 8/80	●Fluorinated Ethylene Propylene (FEP); ●Stripes: Pantone 542 blue; ●X-ray contrast agent: 24% BaSO₄ (pre-mixed in main material)
C22937S	Membrane	●LSR (Liquid silicone rubber)
C31757	PTFE (Polytetrafluoroethylene) Etching Agent	●Etching Agent for PTFE
J1040081		
C21834	PEG handle piece Luer pos.	●PBT (Polybutylenterephthalate) ●Resin ester polymer (carrier material) + pigment dispersion
C21835	PEG LL pos	●PBT (Polybutylenterephthalate) ●Resin ester polymer (carrier material) + pigment dispersion
C22125	LL neg. (cap)	●PBT (Polybutylenterephthalate) ●Resin ester polymer (carrier material) + pigment dispersion
C22126	Lug for PEG connection (cap)	●PA 6 (Polyamid) ●Resin ester polymer (carrier material) + pigment dispersion
C22127	Handle piece (cap)	●PBT (Polybutylenterephthalate) ●Resin ester polymer (carrier material) + pigment dispersion
J1040088		
C22201	Gastric retention disc	●PDMS (Polydimethylesiloxane) (Silicone elastomer) containing vinyl groups, with Platinum catalyst and inhibitors ●Bariumsulfate ●PDMS (Polydimethylsiloxane) containing Vinyl-groups
C22503	Ethanol absolute	●Ethanol absolute
C22559	Glue	●ECA (Ethyl-2-Cyanoacrylate); Hydrochinon

Article code	Article name	Material
J1040063	tube injection CH 15	●TPU (thermoplastic Polyurethane), aromatic-based ●PUR (Polyurethane) ●tube printing ink: Ketone based colourant (cyclohexanone, butanone)
J1040090		
C21891	External fixation plate	●PDMS (Polydimethylsiloxane) (Silicone elastomer) containing vinyl groups, with Platinum catalyst and inhibitors ●Bariumsulphate ● PDMS (Polydimethylsiloxane) containing vinyl groups
C21893	Fixation clamp	●POM (Polyoxymethylene) or Acetal-Copolymer, base resin, additives ● Resin ester polymer (carrier material) + pigment dispersion
J1040104		
C21872	Proximal cone	●TPU (thermoplastic Polyurethane), Polyester-based ●Carrier material: resin-ester polymer with pigment additive
C22559	Glue	●ECA (Ethyl-2-Cyanoacrylate); Hydrochinon
C22691	Fixation loop, OD 0.4 ± 0.01mm	●Core: PA 6.6 (Polyamide 6.6) ●Coating: PA 6 (Polyamide 6)
J1040207		
C21857	Fixation aid PEG	●PBT (Polybutylenterephthalate) ●Resin ester polymer (carrier material) + pigment dispersion
C21859	Fixation ring	●TPE (Thermoplastic Elastomer) ● color: white
C21861	Fixation screw	●PBT (Polybutyleneterephthalate) ●Resin ester polymer (carrier material) + pigment dispersion
Packaging		
	Label Sticker 50 x 22 mm	●Paper and Adhesive
	Adhesive label white 150 x 80 mm	● Upper basis paper: white Glassine paper B700-473, 55µm, 62 g/m² ●Writing paper: satined article 1080, 74 g/m², 71µm The writing paper is coated from both sides with the adhesive ●Adhesive: Permanent acrylic dispersion, 958 ●Bottom basis paper: Glass
	Label sticker 140 x 175 mm	●Paper and Adhesive (acrylate dispersion)
	Natron Paper 1000 x 750 mm	●Color brown

Article code	Article name	Material
	Imprinted envelope	● Printed paper (Velin 80 white, printed colour P 294)
	Instructions for use	● Paper
	Stretch film	● LLDPE (linear low density polyethylene)
		● LDPE (low density polyethylene)
	Cardboard 290 x 187 x 205 mm	● Corrugated cardboard
	Compound foil 30/70 µm	● PA (Polyamide)/PE (Polyethylene)
	Automat reel without print	● Paper
Sterilization		
	Indicator tape for EO sterilization	● Process indicator tape for Ethylene Oxide (EO) sterilization
	Sterilization gas	● Ethylene oxide