

Article code	7901191AP
Article name acc. DoC	Freka Intestinal Tube FR 9 for PEG FR 15
Class acc. Directive 93/42/EEC	IIb

Article code	Article name	Material
C11515	Solvent	●MEK (Methylethylketone)
C11553	Solvent	●CHN (Cyclohexanone)
C21867	Closure cap Luer Lock pos.	●TPU (thermoplastic Polyurethane), Polyether-based ●Preparation of colour agent in a universal carrier
C21898	End piece (flexible tip)	●PUR (Polyurethane)
C22537	Guide wire, 1.14 x 1460 mm	●Stainless steel, coating: PTFE (Polytetrafluoroethylene)
C22563	Mini tube 2.0 x 0.18 x 11.5 mm	●Stainless steel 1.4541
C22572	Lubricant	●PDMS (Polydimethylsiloxane)
C22918	locking screw FR 15	●PBT (Polybutyleneterephthalate), color batch: universal Masterbatch ●PBT (Polybutyleneterephthalate) ●TPE (thermoplastic elastomere)
C30654	Luer Lock pos.	●TPU (thermoplastic Polyurethane), Polyester-based ●TPU (thermoplastic Polyurethane), Polyester-based ●Preparation of colour agent in a universal carrier
C30969	Tube 3.5 x 0.58 mm	●PUR (Polyurethane) and PUR (Polyurethane) with 60% BaSO4 + Colourant (white)
J1040085		
C21835	PEG LL pos	●PBT (Polybutylenterephthalate) ●Resin ester polymer (carrier material) + pigment dispersion
C21933	PEG handle piece Luer pos.	●PBT (Polybutylenterephthalate) ●Resin ester polymer (carrier material) + pigment dispersion
C22125	LL neg. (cap)	●PBT (Polybutylenterephthalate) ●Resin ester polymer (carrier material) + pigment dispersion
C22126	Lug for PEG connection (cap)	●PA 6 (Polyamid) ●Resin ester polymer (carrier material) + pigment dispersion
C22127	Handle piece (cap)	●PBT (Polybutylenterephthalate) ●Resin ester polymer (carrier material) + pigment dispersion
C22568	Mini tube 2.2 x 0.18 x 11.5 mm	●Stainless steel
C22596	Glue	●Epoxy Resin (Reaction mass of trientine and trientine, mono- and di-propoxylated)

Article code	Article name	Material
C22597	Glue hardener	● Modified Polyamine (2,2'-[(1-methylethylidene) bis(4,1-phenyleneoxymethylene)] bisoxirane , diethyl phthalate)
J1040097		
C21835	PEG LL pos	● PBT (Polybutylenterephthalate) ● Resin ester polymer (carrier material) + pigment dispersion
C22010	Y-piece	● PBT (Polybutylenterephthalate) ● Resin ester polymer (carrier material) + pigment dispersion
C22125	LL neg. (cap)	● PBT (Polybutylenterephthalate) ● Resin ester polymer (carrier material) + pigment dispersion
C22126	Lug for PEG connection (cap)	● PA 6 (Polyamid) ● Resin ester polymer (carrier material) + pigment dispersion
C22127	Handle piece (cap)	● PBT (Polybutylenterephthalate) ● Resin ester polymer (carrier material) + pigment dispersion
C22559	Glue	● ECA (Ethyl-2-Cyanoacrylate); Hydrochinon
C30660	PEG LL pos.	● PBT (Polybutylenterephthalate) ● Resin ester polymer (carrier material) + pigment additive
J1040098		
C21846	Membrane	● TPE (Thermoplastic Elastomer) ● white
C21847	Fixation screw	● PBT (Polybutylenterephthalate) ● Resin ester polymer (carrier material) + pigment dispersion
C21848	Adaptor interior part	● PBT (Polybutylenterephthalate) ● Carrier: resin-ester polymer with pigment dispersion
C21849	Adaptor exterior part	● POM (Polyoxymethylene) or Acetal-Copolymer, base resin, additives ● Carrier: resin-ester polymer with pigment dispersion
C21858	Fixation ring	● TPE (Thermoplastic Elastomer) ● color: white
C22559	Glue	● ECA (Ethyl-2-Cyanoacrylate); Hydrochinon
J1040207		
C21857	Fixation aid PEG	● PBT (Polybutylenterephthalate) ● Resin ester polymer (carrier material) + pigment dispersion
C21859	Fixation ring	● TPE (Thermoplastic Elastomer) ● color: white

Article code	Article name	Material
C21861	Fixation screw	●PBT (Polybutyleneterephthalate) ●Resin ester polymer (carrier material) + pigment dispersion
J1060207		
C30979	Tube 1.9 x 0.51	●PUR (Polyurethane) and PUR + 60% BaSO4 (one x-ray contrast strip)
C31534	Tube printing ink	●Butanone ●Solvent Black 29 ●Ethanol ●Cellulose Nitrate ●Isopropyle Acetate ●2-Naphtole
C31535	Tube printing ink thinner	●Butanone ●Acetone
J1060209		
C21988	Adaptor with Luer Lock neg.	●PA 6 (Polyamide)
C21989	Seal ring 3.1 x 6.5 x 1.6 mm	●PE (Polyethylene)
C21990	Adaptor insertion	●SEBS (Styrene-ethylene/ butylene-styrene)
C21991	Y-piece	●PC (Polycarbonate)
Packaging		
	Compound foil 70 µm	●PA (Polyamide) / adhesive layer / PE (Polyethylene)
	Peel paper	●Paper 60 g/m2
	Label sticker 105 x 144 mm	●Paper
	Label Sticker 50 x 22 mm	●Paper and Adhesive
	Adhesive label white 150 x 80 mm	● Upper basis paper: white Glassine paper B700-473, 55µm, 62 g/m² ●Writing paper: satined article 1080, 74 g/m², 71µm The writing paper is coated from both sides with the adhesive ●Adhesive: Permanent acrylic dispersion, 958 ●Bottom basis paper: Glass
	Instructions for use	●Printed paper
	Self adhesive tape	●PVC (Polyvinylchloride) as carrier foil with natural rubber as adhesive
	Natron Paper 1000 x 750 mm	●Color brown
	Stretch film	●LLDPE (linear low density polyethylene)
	Stretch film foil plate	●LDPE (low density polyethylene)
	Folding box 782 x 595 x 360 mm	●Corrugated cardboard
	Cardboard 253 x 150 x 85 mm	●Corrugated cardboard, color white
	Blue ink	●Printing ink
Sterilization		

Article code	Article name	Material
	Indicator tape for EO sterilization	● Process indicator tape for Ethylene Oxide (EO) sterilization
	Sterilization gas	● Ethylene oxide