

Article code	7980111
Article name acc. DoC	Freka® Tube FR 15, 100 cm, F
Class acc. Directive 93/42/EEC	IIa

Article code	Article name	Material
C10575	Mandrin plug	●PC (Polycarbonate) ●Resin ester polymer (carrier material)
C11515	Solvent	●MEK (Methylethylketone)
C11553	Solvent	●CHN (Cyclohexanone)
C21986	Tube tip	●PUR (Polyurethane) with 10% BaSO₄ ●Lubricant (Amid wax)
C22094	Funnel adaptor for mandrin	●MABS (Methymethacrylate Acrylonitrile Butadiene Styrene) ●Resin ester polymer (carrier material)
C22501	Aceton - for analysis	●Acetone (2-Propanone, Dimethyl ketone)
C22559	Glue	●ECA (Ethyl-2-Cyanoacrylate); Hydrochinon (0.01- < 0.1% w/w)
C22572	Lubricant	●PDMS (Polydimethylsiloxane)
C22631	Mandrin 1.7 x 1002 mm	●Stainless steel
C22704	Fixation tape	●polyester fleece non-woven ●polyacrylate adhesive ●siliconized paper liner
C22914	Mandrin 1.7 x 1002 mm, white	●C22631 Sondenmandrin/Spiral L=1002 mm rd 1.7 ●C22094 Funnel mandrin attachment ●C10575 Steel mandrin end stopper
C30036	Connecting funnel	●TPU (thermoplastic Polyurethane) ●Resin ester polymer (carrier material)
C30969	Tube 3.5 x 0.58 mm	●PUR (Polyurethane) and PUR (Polyurethane) with 60% BaSO₄ + Colourant (white)
C31534	Tube printing ink	●Butanone ●Solvent Black ●Ethanol ●Cellulose Nitrate ●Isopropyle Acetate ●2-Naphtole
C31535	Tube printing ink thinner	●Butanone ●Acetone
C60968	TPU-tube 3,50 x, 0,50 mm with 1 RKS	Tube: ●TPU (thermoplastic Polyurethane), Polyester-based; color: natural; X-ray contrast stripe: ●TPU (thermoplastic Polyurethane), Polyester-based + 60% BaSO₄; color: white;

Article code	Article name	Material
C90183	Elastollan 1185 A 10 000	●TPU (thermoplastic Polyurethane), Polyether-based
Packaging		
	Edge protection	●corrugated cardboard
	Compound foil 70 µm	●PA (Polyamide) / adhesive layer / PE (Polyethylene)
	Label sticker 145 x 94 mm	●80 g/sqm adhesive paper (adhesive based on acrylic)
	Natron Paper 1000 x 750 mm	●Color brown
	Stretch film	●LLDPE (linear low density polyethylene)
	Stretch film foil plate	●LDPE (low density polyethylene)
	Cardboard 249 x 148 x 100 mm	●Corrugated cardboard, color brown inside and outside
	Cardboard 249 x 148 x 100 mm	●Corrugated cardboard, color brown ●Printing FK Logo: Blue- Pantone 300U
	Black ink	●Printing ink
	Blue ink	●Printing ink
	Packaging medical paper 60 g/m2	●Paper
Sterilization		
	Indicator tape for EO sterilization	●Process indicator tape for Ethylene Oxide (EO) sterilization
	Sterilization gas	●Ethylene oxide