

Material List
Finished Product: 7980321



Internal code	Designation	Polymer
C22540	Guide wire 0,035x 3000 mm	●Dispenser-Tube ●Mandrin ●Insertion guide: LDPE (low density polyethylene)
<u>Packaging</u>		
C11008	Folding box 394 x 288 x 324 mm	●Corrugated cardboard, color brown
C11517	Adhesive tape paper	●gummed paper tapes Polyeste filament yarn for strengthening
C12055	Tamper-evident safety seal	●PP (Polypropylene) -Film Transparent- 60 µm Top Coat ●Adhesive: Permanent -958 ●Base: Glassine paper B 700-476
C12182	Roll Label 105 x 148 mm	●Paper (blank white) and Adhesive
C31045	Stretch film	●LLDPE (linear low density polyethylene)
C31046	Stretch film foil plate	●LDPE (low density polyethylene)
C31351	Compound foil 30/70 µm	●PA (Polyamide)/PE (Polyethylene)
C31525	Printing ink, blue	●Mixture (colourant and thinners)
C31526	Thinner for printing ink	●1-Methoxy-2-Propanol (2-Metholymethylethoxy) propanol
C51027	Folded carton box, brown, 260x250x35 mm	● Material: Testliner ● Wave: E
C61300	Automat reel without print	●Paper
<u>Sterilization</u>		
C11525	Indicator tape for EO sterilization	●Process indicator tape for Ethylene Oxide (EO) sterilization
C31500	Sterilization gas, Poland	●100% Ethylene oxide, liquefied grade one (GAZ 4-134) ●alternative: > 99.9% Ethylene oxide, CAS no. 75-21-8