

Article code	7980961
Article name acc. DoC	Freka® Paed Tube FR 6.5, 50 cm, F
Class acc. Directive 93/42/EEC	IIa

Article code	Article name	Material
C11515	Solvent	●MEK (Methylethylketone)
C11553	Solvent	●CHN (Cyclohexanone)
C21874	Connecting funnel	●TPU (thermoplastic Polyurethane)
C22572	Lubricant	●PDMS (Polydimethylsiloxane)
C30982	Tube 1.55 x 2.38 mm	●PUR (Polyurethane) and PUR (Polyurethane) with 60% BaSO4 + Colourant (white)
C31534	Tube printing ink	●Butanone ●Solvent Black ●Ethanol ●Cellulose Nitrate ●Isopropyle Acetate ●2-Naphtole
C31535	Tube printing ink thinner	●Butanone ●Acetone
Packaging		
	Compound foil 70 µm	●PA (Polyamide) / adhesive layer / PE (Polyethylene)
	Peel paper	●Paper 60 g/m2
	Self-adhesive tape, white	●PP (Polypropylene) ●adhesive: natural rubber
	Natron Paper 1000 x 750 mm	●Color brown
	Carton sticker	●80 g/qm adhesive paper
	Instruction leaflet 120 x 205 mm	●50 g/qm white-frost paper
	Stretch film	●LLDPE (linear low density polyethylene)
	Stretch film foil plate	●LDPE (low density polyethylene)
	Corrugated folding box 380 x 260 x 160 mm	●Corrugated cardboard
	Cardboard 253 x 150 x 85 mm	●Corrugated cardboard, color white
	Protective Layer 200 x 370 mm	●cardboard
	Black ink	●Printing ink
Sterilization		
	Indicator tape for EO sterilization	●Process indicator tape for Ethylene Oxide (EO) sterilization
	Sterilization gas	●Ethylene oxide