

Article code	7980961
Article name acc. DoC	Freka® Paed Tube FR 6.5, 50 cm, F
Class acc. Directive 93/42/EEC	IIa

Article code	Article name	Material
C11515	Solvent	●MEK (Methylethylketone)
C11553	Solvent	●CHN (Cyclohexanone)
C21874	Connecting funnel	●C90183 TPU (thermoplastic Polyurethane) ●Lubricant: Amide wax
C22572	Lubricant	●PDMS (Polydimethylsiloxane)
C30982	Tube 1.55 x 2.38 mm	●PUR (Polyurethane) and PUR (Polyurethane) with 60% BaSO4 + Colourant (white)
C31534	Tube printing ink	●Butanone ●Solvent Black 29 ●Ethanol ●Cellulose Nitrate ●Isopropyle Acetate ●2-Naphtole
C31535	Tube printing ink thinner	●Butanone ●Acetone
Packaging		
	Compound foil 70 µm	●PA (Polyamide) / adhesive layer / PE (Polyethylene)
C11302	Peel paper	●Paper 60 g/m2
C22621	Self-adhesive tape, white	●PP (Polypropylene) ●adhesive: natural rubber
C22645	Natron Paper 1000 x 750 mm	●Color brown
C22890	Carton sticker	●80 g/qm adhesive paper
	Instruction leaflet 120 x 205 mm	●50 g/qm white-frost paper
C31045	Stretch film	●LLDPE (linear low density polyethylene)
C31046	Stretch film foil plate	●LDPE (low density polyethylene)
C31062	Corrugated folding box 380 x 260 x 160 mm	●Corrugated cardboard
C31063	Cardboard 253 x 150 x 85 mm	●Corrugated cardboard, color white
C31077	Protective Layer 200 x 370 mm	●cardboard
	Black ink	●Printing ink
Sterilization		
C11525	Indicator tape for EO sterilization	●Process indicator tape for Ethylene Oxide (EO) sterilization

Article code	Article name	Material
C31500	Sterilization gas, Poland	●100% Ethylene oxide, liquefied grade one (GAZ 4-134) ●alternative: > 99.9% Ethylene oxide