

Article code	7981367
Article name acc. DoC	Insert for Freka® GastroTube FR 15, ENFit
Class acc. MDR 2017/745	I unsterile

Article code	Article name	Material
C30419	ENFit male insert for GastroTube	●MABS (Methylmethacrylate-Acrylonitrile-Butadiene-Styrene), ●Masterbatch violett (RAL4005)
Packaging		
	Roll ENFit label 105×148.8 mm	●Paper (blank white with violette bar) and Adhesive
	Additional information	●Paper
	Natron Paper 1000 x 750 mm	●Color brown
	Stretch film	●LLDPE (linear low density polyethylene)
	Stretch film foil plate	●LDPE (low density polyethylene)
	Cardboard 290 x 187 x 205 mm	●Corrugated cardboard
	Zip lock bag	●PE (Polyethylene)