

Article code	7981388
Article name acc. DoC	Freka® Repair Set for Freka® FCJ, ENFit
Class acc. MDR 2017/745	I

Article code	Article name	Material
C21941	Tube clamp with centering, white	●POM (Polyoxymethylene) ●Color Carrier: resin-ester polymer with pigment dispersion
J1040198		
C22563	Mini tube 2.0 x 0.18 x 11.5 mm	●Stainless steel 1.4541
C22596	Glue	●Epoxy Resin (Reaction mass of trientine and trientine, mono- and di-propoxylated EC no. 942-835-1)
C22597	Glue hardener	●Modified Polyamine (2,2'-[(1-methylethyldiene) bis(4,1-phenyleneoxymethylene)] bisoxirane , diethyl phthalate)
C30823	ENFit female cap	●HDPE (high density Polyethylene) ●colourant violet
C30894	ENFit male - PEG connector	●PBT (Polybutylenterephthalate) ●Resin ester polymer (carrier material)
C30897	Strap for ENFit male connector	●LDPE (Low density polyethylene)
J1040303		
C21857	Fixation aid PEG	●PBT (Polybutylenterephthalate) ●Resin ester polymer (carrier material) + pigment dispersion
C21858	Fixation ring	●TPE (Thermoplastic Elastomer) ● color: white
C21868	Fixation screw	●PBT (Polybutyleneterephthalate) ●Resin ester polymer (carrier material) + pigment dispersion
Packaging		
	Roll ENFit label 105×148.8 mm	●Paper (blank white with violette bar) and Adhesive
	IFU Freka Repair Set for Freka FCJ	●Paper
	Natron Paper 1000 x 750 mm	●Color brown
	Stretch film	●LLDPE (linear low density polyethylene)
	Stretch film foil plate	●LDPE (low density polyethylene)
	Cardboard 290 x 187 x 205 mm	●Corrugated cardboard
	Zip lock bag	●PE (Polyethylene)