

Article code	7981396
Article name acc. DoC	Freka® Y-Connector for Freka® PEG FR 20/12, ENFit
Class acc. MDR 2017/745	I unsterile

Article code	Article name	Material
J1040135		
C11503	Solvent	●THF (Tetrahydrofuran), also 1,4-Epoxy butane
C11553	Solvent	●CHN (Cyclohexanone)
C21811	Connecting funnel	●TPU (thermoplastic Polyurethane), Polyether-based ●Carrier: resin-ester polymer with pigment dispersion ●Lubricant (Amid wax)
C21833	Sealing body for click adaptor	●PBT (Polybutylenterephthalate) ●Resin ester polymer (carrier material) + pigment dispersion
C30808	ENLock Y-piece	●PBT (Polybutylenterephthalate) ●Resin ester polymer (carrier material) + pigment dispersion
J1040300		
C21839	Fixation ring Ø6,5 mm CH20	●TPE (Thermoplastic Elastomer) ●Lubricant: Amide wax
C21857	Fixation aid PEG	●PBT (Polybutylenterephthalate) ●Resin ester polymer (carrier material) + pigment dispersion
C30812	Fixation screw	●PBT (Polybutylenterephthalate) ●Resin ester polymer (carrier material) + pigment dispersion
Packaging		
	Roll ENFit label 105×148.8 mm	●Paper (blank white with violette bar) and Adhesive
	Additional information	Paper
	Natron Paper 1000 x 750 mm	●Color brown
	Stretch film	●LLDPE (linear low density polyethylene)
	Stretch film foil plate	●LDPE (low density polyethylene)
	Cardboard 290 x 187 x 205 mm	●Corrugated cardboard
	Zip lock bag	●PE (Polyethylene)