

Article code	7981823AP
Article name acc. DoC	Freka® Sonde CH 12, 120 cm, ENLock
Class acc. Directive 93/42/EEC	IIa

Article code	Article name	Material
C11503	Solvent	●THF (Tetrahydrofuran), also 1,4-Epoxy butane
C11553	Solvent	●CHN (Cyclohexanone)
C22501	Aceton - for analysis	●Acetone (2-Propanone, Dimethyl ketone)
C22572	Lubricant	●PDMS (Polydimethylsiloxane)
C22704	Fixation tape	●polyester fleece non-woven ●polyacrylate adhesive ●siliconized paper liner
C30684	Connecting funnel ENLock	●TPU (thermoplastic Polyurethane) ●Resin ester polymer (carrier material) + pigment dispersion ●TPU (thermoplastic Polyurethane)
C30686	Cap	●POE (Polyolefine elastomers) ●lubricant copolymer of ethylene + octene ●(carrier material) + pigment dispersion
C30691	ENLock female Y	●TPU (thermoplastic Polyurethane) ●Resin ester polymer (carrier material) + pigment dispersion
C30729	Tube Tip FR 12	●PUR (Polyurethane) with 10% BaSO4 ●Lubricant (Amidwachs) ●Slip additive Colorslip 103 100 % carboxylic acid ester
C30913	Tube 2.9 x 4 mm	●PUR (Polyurethane) + PUR with 60% BaSO4 ●Colourant (white)
C31534	Tube printing ink	●Butanone ●Solvent Black 29 ●Ethanol ●Cellulose Nitrate ●Isopropyle Acetate ●2-Naphtole
C31535	Tube printing ink thinner	●Butanone ●Acetone
J1040122		
C22629	Mandrin 1.7 x 1225 mm	●Stainless steel
C30815	Handle Piece	● MABS (Methylmethacrylate Acrylnitrile Butadiene Styrene) ●Resin ester polymer (carrier material) + pigment dispersion
C30819	Glue	●superglue based on Alkoxyethyl cyanoacrylate

Article code	Article name	Material
Packaging		
	Edge protection	●corrugated cardboard
	Compound foil 70 µm	●PA (Polyamide) / adhesive layer / PE (Polyethylene)
	White paper for IFU printing	●100% wood pulp
	Label sticker 145 x 94 mm	●80 g/sqm adhesive paper (adhesive based on acrylic)
	Natron Paper 1000 x 750 mm	●Color brown
	Stretch film	●LLDPE (linear low density polyethylene)
	Stretch film foil plate	●LDPE (low density polyethylene)
	Cardboard 249 x 148 x 100 mm	●Corrugated cardboard, color brown ●Printing FK Logo: Blue- Pantone 300U
	Black ink	●Printing ink
	Packaging medical paper 60 g/m2	●Paper
Sterilization		
	Indicator tape for EO sterilization	●Process indicator tape for Ethylene Oxide (EO) sterilization
	Sterilization gas	●Ethylene oxide