

Article code	7981826
Article name acc. DoC	Freka Tube FR 8, 80cm, ENFit
Class acc. Directive 93/42/EEC	IIa

Article code	Article name	Material
C11515	Solvent	●MEK (Methylethylketone)
C11553	Solvent	●CHN (Cyclohexanone)
C21948	Tube tip	●UR (Polyurethane) with 10% BaSO ₄
C22501	Aceton - for analysis	●Acetone (2-Propanone, Dimethyl ketone)
C22572	Lubricant	●PDMS (Polydimethylsiloxane)
C22704	Fixation tape	●polyester fleece non-woven ●polyacrylate adhesive ●siliconized paper liner
C30352	ProNeo Soft Cap, white	●POE (Polyolefine elastomers) ●lubricant copolymer of ethylene + octene ●colourant - white
C30777	ENFit Insert for Y port catheter connector	● MABS (Methylmethacrylate Acrylnitrile Butadiene Styrene) ●blue lilac MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene-basis + pigment)
	ENFit Y port 8 FR	●TPU (thermoplastic Polyurethane) ●mould release agent (2%): TPU-(thermoplastic Polyurethane)-based anti-blocking agent ●(former E6544) colourant violett (2%) masterbatch, Carrier: thermoplastic polyurethane with inorganic and
C30981	Tube 1.9 x 2.85 mm	●PUR (Polyurethane) and PUR (Polyurethane) with 60% BaSO ₄ + Colourant (white)
C31534	Tube printing ink	●Butanone ●Solvent Black 29 ●Ethanol ●Cellulose Nitrate ●Isopropyle Acetate ●2-Naphtole
C31535	Tube printing ink thinner	●Butanone ●Acetone
J1040259		
	Mandrin 1.37 x 805 mm	●Stainless steel
C30819	Glue	●superglue based on Alkoxyethyl cyanoacrylate

Article code	Article name	Material
C30885	Guide wire grip ENFit, 1.3 mm	●MABS (Methylmethacrylate Acrylnitrile Butadiene Styrene) ● Resin ester polymer (carrier material), colourant white
C31786	Mandrin L = 805 mm	●Stainless steel 1.4310
Packaging		
C11096	Edge protection	●corrugated cardboard
	Compound foil 70 µm	●PA (Polyamide) / adhesive layer / PE (Polyethylene)
C11302	Peel paper	●Paper 60 g/m²
C12120	Label sticker 145 x 94 mm	FKB: ●80 g/sqm adhesive paper (adhesive based on acrylic) FKN: ●high gloss paper / adhesive / liner
C22645	Natron Paper 1000 x 750 mm	●Color brown
C31045	Stretch film	●LLDPE (linear low density polyethylene)
C31046	Stretch film foil plate	●LDPE (low density polyethylene)
	Cardboard 249 x 148 x 100 mm	●Corrugated cardboard, color white
C31561	Black ink	●Printing ink
C31562	Blue ink	●Printing ink
Sterilization		
	Indicator tape for EO sterilization	●Process indicator tape for Ethylene Oxide (EO) sterilization
C31500	Sterilization gas, Poland	●100% Ethylene oxide, liquefied grade one (GAZ 4-134) ●alternative: > 99.9% Ethylene oxide