

Article code	7981827
Article name acc. DoC	Freka Tube FR 8, 120cm, ENFit
Class acc. Directive 93/42/EEC	IIa

Article code	Article name	Material
C11515	Solvent	•MEK (Methylethylketone)
C11553	Solvent	•CHN (Cyclohexanone)
C21948	Tube tip	•PUR (Polyurethane) with 10% BaSO ₄
C22501	Aceton - for analysis	•Acetone (2-Propanone, Dimethyl ketone)
C22572	Lubricant	•PDMS (Polydimethylsiloxane)
C22704	Fixation tape	•polyester fleece non-woven •polyacrylate adhesive •siliconized paper liner
C30352	ProNeo Soft Cap, white	•POE (Polyolefine elastomers) •lubricant copolymer of ethylene + octene •colourant - white
C30777	ENFit Insert for Y port catheter connector	• MABS (Methylmethacrylate Acrylnitrile Butadiene Styrene) •blue lilac MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene-basis + pigment)
	ENFit Y port 8 FR	•TPU (thermoplastic Polyurethane) •mould release agent (2%): TPU- (thermoplastic Polyurethane)-based anti-blocking agent •(former E6544) colourant violett (2%) masterbatch, Carrier: thermoplastic polyurethane with inorganic and
C30981	Tube 1.9 x 2.85 mm	•PUR (Polyurethane) and PUR (Polyurethane) with 60% BaSO ₄ + Colourant (white)
C31534	Tube printing ink	•Butanone •Solvent Black 29 •Ethanol •Cellulose Nitrate •Isopropyle Acetate •2-Naphtole
C31535	Tube printing ink thinner	•Butanone •Acetone
J1040260		
	Mandrin 1.37 x 1205 mm	•Stainless steel
C30819	Glue	•superglue based on Alkoxyethyl cyanoacrylate

Article code	Article name	Material
C30885	Guide wire grip ENFit, 1.3 mm	• MABS (Methylmethacrylate Acrylnitrile Butadiene Styrene) • Resin ester polymer (carrier material), colourant white
C31789	Mandrin L = 1205 mm	• Stainless steel 1.4310
Packaging		
C11096	Edge protection	• corrugated cardboard
	Compound foil 70 µm	• PA (Polyamide) / adhesive layer / PE (Polyethylene)
C11302	Peel paper	• Paper 60 g/m²
C12120	Label sticker 145 x 94 mm	FKB: • 80 g/sqm adhesive paper (adhesive based on acrylic) FKN: • high gloss paper / adhesive / liner
C22645	Natron Paper 1000 x 750 mm	• Color brown
C31045	Stretch film	• LLDPE (linear low density polyethylene)
C31046	Stretch film foil plate	• LDPE (low density polyethylene)
	Cardboard 249 x 148 x 100 mm	• Corrugated cardboard, color white
C31561	Black ink	• Printing ink
C31562	Blue ink	• Printing ink
Sterilization		
	Indicator tape for EO sterilization	• Process indicator tape for Ethylene Oxide (EO) sterilization
C31500	Sterilization gas, Poland	• 100% Ethylene oxide, liquefied grade one (GAZ 4-134) • alternative: > 99.9% Ethylene oxide