

Article code	7981830
Article name acc. DoC	Freka Tube FR 15, 120cm, ENFit
Class acc. Directive 93/42/EEC	IIa

Article code	Article name	Material
C10575	Mandrin plug	●PC (Polycarbonate) ●Resin ester polymer (carrier material)
C11515	Solvent	●MEK (Methylethylketone)
C11553	Solvent	●CHN (Cyclohexanone)
C21986	Tube tip	●PUR (Polyurethane) with 10% BaSO₄ ●Lubricant (Amid wax)
C22501	Aceton - for analysis	●Acetone (2-Propanone, Dimethyl ketone)
C22559	Glue	●ECA (Ethyl-2-Cyanoacrylate); Hydrochinon (0.01- < 0.1% w/w)
C22572	Lubricant	●PDMS (Polydimethylsiloxane)
C22704	Fixation tape	●polyester fleece non-woven ●polyacrylate adhesive ●siliconized paper liner
C30352	ProNeo Soft Cap, white	●POE (Polyolefine elastomers) ●lubricant copolymer of ethylene + octene ●colourant - white
C30777	ENFit Insert for Y port catheter connector	●MABS (Methylmethacrylate Acrylnitrile Butadiene Styrene) ●blue lilac MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene-basis + pigment)
C30781	ENFit Y port 15 FR	●TPU (thermoplastic Polyurethane) ●mould release agent (2%): TPU-(thermoplastic Polyurethane)-based anti-blocking agent ●colourant violett (2%) masterbatch, Carrier: thermoplastic polyurethane with inorganic an
C30889	Guide wire ENFit, 1.7 mm, 1225 mm	●C30888: MABS (Methylmethacrylate Acrylnitrile Butadiene Styrene) + Masterbatch ●C22629 Stainless Steel ●C30819 glue Loctite 4031: Cyanoacrylate
C30969	Tube 3.5 x 0.58 mm	●PUR (Polyurethane) and PUR (Polyurethane) with 60% BaSO₄ + Colourant (white)
C31534	Tube printing ink	●Butanone ●Solvent Black ●Ethanol ●Cellulose Nitrate ●Isopropyle Acetate ●2-Naphtole

Article code	Article name	Material
C31535	Tube printing ink thinner	● Butanone ● Acetone
C60968	TPU-tube 3,50 x, 0,50 mm with 1 RKS	Tube: ● TPU (thermoplastic Polyurethane), Polyester-based; color: natural; X-ray contrast stripe: ● TPU (thermoplastic Polyurethane), Polyester-based + 60% BaSO₄; color: white;
Packaging		
	Edge protection	● corrugated cardboard
	Compound foil 70 µm	● PA (Polyamide) / adhesive layer / PE (Polyethylene)
	Peel paper	● Paper 60 g/m²
	Label sticker 145 x 94 mm	● 80 g/sqm adhesive paper (adhesive based on acrylic)
	Natron Paper 1000 x 750 mm	● Color brown
	Stretch film	● LLDPE (linear low density polyethylene)
	Stretch film foil plate	● LDPE (low density polyethylene)
	Cardboard 249 x 148 x 100 mm	● Corrugated cardboard, color white
	Black ink	● Printing ink
	Blue ink	● Printing ink
Sterilization		
	Indicator tape for EO sterilization	● Process indicator tape for Ethylene Oxide (EO) sterilization
	Sterilization gas	● Ethylene oxide