

Article code	7981832
Article name acc. DoC	Freka® Endolumina, ENFit
Class acc. Directive 93/42/EEC	IIa

Article code	Article name	Material
9000784	Freka Connect ENFit / ProNeo syringe	●Syringe Barrel: PP (Polypropylene) ●Syringe Plunger: PP (Polypropylene) ●Syringe Gasket - Isoprene Rubber ●Lubricant - medical silicone oil ●Print - PPE (Polyphenylenether) Ink Series
C21857	Fixation aid PEG	●PBT (Polybutylenterephthalate) ●Resin ester polymer (carrier material) + pigment dispersion
C21858	Fixation ring	●TPE (Thermoplastic Elastomer) ● color: white
C21869	Fixation screw	●PBT (Polybutylenterephthalate) ●Resin ester polymer (carrier material)
C22572	Lubricant	●PDMS (Polydimethylsiloxane)
C22704	Fixation tape	●polyester fleece non-woven ●polyacrylate adhesive ●siliconized paper liner
C22907	Assembly screwing aid CH8, grey	●C21858: TPE (thermoplastic elastomere) ●C21869: PBT (Polybutyleneterephthalate) ●C21857: PBT (Polybutyleneterephthalate) + Masterbatch
C30765	Endo-Schale 200x200 mm, ENFit	●PETG (Polyethyleneterephthalate Glycol)
C30834	ENFit female oral / LL	●POM (Polyoxymethylene)
C31534	Tube printing ink	●Butanone ●Solvent Black 29 ●Ethanol ●Cellulose Nitrate ●Isopropyle Acetate ●2-Naphtole
C31535	Tube printing ink thinner	●Butanone ●Acetone
C60911	3-layer tube 1.8 x 0.4 x 2.6 mm	●inner layer: Thermoplastic Polyurethane (TPU) Compound + slip additive ●middle layer: Thermoplastic Polyurethane (TPU) with barium sulfate ●outer layer: Thermoplastic Polyurethane (TPU) Compound + slip additive
J1040027		

Article code	Article name	Material
C22174	MCT-oil and vials	●Oil: medium-chain triglycerides - caproic acid, caprylic acid, capric acid, lauric acid ●Vials: LDPE (low density polyethylene)
C22636	Seal bags 145 x 60 mm	●PES (Polyester) ●ALU (Aluminium) foil ●PE (Polyethylene)
J1040108		
C11515	Solvent	●MEK (Methylethylketone)
C11553	Solvent	●CHN (Cyclohexanone)
C21879	Connecting funnel	●TPU (thermoplastic Polyurethane) ●Preparation of a colour agent in a universal carrier ●Lubricant (Amid wax) ●TPU-(thermoplastic Polyurethane)-based anti-blocking agent
C22503	Ethanol absolute	●Ethanol absolute
C22572	Lubricant	●PDMS (Polydimethylsiloxane)
C30975	Tube 3.25 x 0.9 x 350 mm	●PVC (Polyvinylchloride), plasticizer DEHP, Colourant
J1040194		
C22564	Mini tube 2 x 0.18 x 11.5 mm	●stainless steel 1.4541
C22596	Glue	●Epoxy Resin (Reaction mass of trientine and trientine, mono- and di-propoxylated)
C22597	Glue hardener	●Modified Polyamine (2,2'-[(1-methylethyldiene) bis(4,1-phenyleneoxymethylene)] bisoxirane, diethyl phthalate)
C30823	ENFit female cap	●HDPE (high density Polyethylene) ●colourant violet
C30890	ENFit male - PEG connector	●PBT (Polybutylenterephthalate) ●Resin ester polymer (carrier material)
C30897	Strap for ENFit male connector	●LDPE (Low density polyethylene)
Packaging		
	Adhesive label white 150 x 80 mm	● Upper basis paper: white Glassine paper, 55µm, 62 g/m² ●Writing paper: satined, 74 g/m², 71µm The writing paper is coated from both sides with the adhesive ●Adhesive: Permanent acrylic dispersion ●Bottom basis paper: Glass
	Instructions leaflet	●Paper
	Self adhesive tape	●PVC (Polyvinylchloride) as carrier foil with natural rubber as adhesive
	Label sticker 160 x 115 mm	●Paper
	Stretch film	●LLDPE (linear low density polyethylene)

Article code	Article name	Material
	Stretch film foil plate	●LDPE (low density polyethylene)
	Folding box 782 x 595 x 360 mm	●Corrugated cardboard
	Carton 325 x 275 x 90 mm	●Cardboard
	Compound foil 30/70 µm	●PA (Polyamide)/PE (Polyethylene)
	Automat reel without print	●Paper
	Imprinted envelope	●Paper
Sterilization		
	Indicator tape for EO sterilization	●Process indicator tape for Ethylene Oxide (EO) sterilization
	Sterilization gas	●Ethylene oxide