

Article code	7981840
Article name acc. DoC	Freka® Tube FR 8 linear, 120 cm, ENFit
Class acc. Directive 93/42/EEC	IIa

Article code	Article name	Material
C11515	Solvent	●MEK (Methylethylketone)
C11553	Solvent	●CHN (Cyclohexanone)
C21948	Tube tip	●PUR (Polyurethane) with 10% BaSO4
C22501	Aceton - for analysis	●Acetone (2-Propanone, Dimethyl ketone)
C22572	Lubricant	●PDMS (Polydimethylsiloxane)
C22704	Fixation tape	●polyester fleece non-woven ●polyacrylate adhesive ●siliconized paper liner
C30823	ENFit female cap	●HDPE (high density Polyethylene) ●colourant violet
C30824	Strap for ENFit cap, ring-Ø 10.5 mm	●LDPE (low density polyethylene), transparent
C30841	ENFit male 8 FR	●PET (Polyethyleneterephthalate) (Copolyester) ●colourant violet
C30981	Tube 1.9 x 2.85 mm	●PUR (Polyurethane) and PUR (Polyurethane) with 60% BaSO4 + Colourant (white)
C31534	Tube printing ink	●Butanone ●Solvent Black ●Ethanol ●Cellulose Nitrate ●Isopropyle Acetate ●2-Naphtole
C31535	Tube printing ink thinner	●Butanone ●Acetone
Packaging		
	Edge protection	●corrugated cardboard
	Compound foil 70 µm	●PA (Polyamide) / adhesive layer / PE (Polyethylene)
	Peel paper	●Paper 60 g/m2
	Label sticker 145 x 94 mm	●80 g/sqm adhesive paper (adhesive based on acrylic)
	Natron Paper 1000 x 750 mm	●Color brown
	Stretch film	●LLDPE (linear low density polyethylene)
	Stretch film foil plate	●LDPE (low density polyethylene)
	Cardboard 249 x 148 x 100 mm	●Corrugated cardboard, color brown ●Printing FK Logo: Blue- Pantone 300U
	Black ink	●Printing ink
	Blue ink	●Printing ink

Sterilization		
	Indicator tape for EO sterilization	●Process indicator tape for Ethylene Oxide (EO) sterilization
	Sterilization gas	●Ethylene oxide