

Article code	7981844 / M908003423
Article name acc. DoC	Freka Tube FR 15 linear, 120 cm, ENFit
Class acc. Directive 93/42/EEC	Ila

Article code	Article name	Material
C10571	Mandrin plug	●PC (Polycarbonate) based on Bisphenole A ●Resin ester polymer (carrier material)
C11515	Solvent	●MEK (Methylethylketone)
C11553	Solvent	●CHN (Cyclohexanone)
C21986	Tube tip	●PUR (Polyurethane) with 10% BaSO₄ ●Lubricant (Amid wax)
C22501	Aceton - for analysis	●Acetone (2-Propanone, Dimethyl ketone)
C22559	Glue	●ECA (Ethyl-2-Cyanoacrylate)
C22572	Lubricant	●PDMS (Polydimethylsiloxane)
C22704	Fixation tape	●polyester fleece non-woven ●polyacrylate adhesive ●siliconized paper liner
C30766	ENFit male, 13 FR	●PET (Polyethyleneterephthalate) (Copolyester) ●colourant violet
C30823	ENFit female cap	●HDPE (high density Polyethylene) ●colourant violet
C30824	Strap for ENFit cap, ring-Ø 10.5 mm	●LDPE (low density polyethylene), transparent
C30969	Tube 3.5 x 0.58 mm	●PUR (Polyurethane) and PUR (Polyurethane) with 60% BaSO₄ + Colourant (white)
C31526	Thinner for printing ink	●1-Methoxy-2-Propanol (2-Metholymethylethoxy) propanol
C31527	Printing ink, black	●Stamping ink
C31534	Tube printing ink	●Butanone 70-80% ●Solvent Black 29 <10% ●Ethanol 5-10% ●Cellulose Nitrate 5-10% ●Isopropyle Acetate 5% ●2-Naphtole <0.2%
C31535	Tube printing ink thinner	●Butanone 96-98% ●Acetone 2-4%
Packaging		
	Edge protection	●corrugated cardboard
	Label sticker 145 x 94 mm	●80 g/sqm adhesive paper (adhesive based on acrylic)
	Natron Paper 1000 x 750 mm	●Color brown
	Stretch film	●LLDPE (linear low density polyethylene)
	Stretch film foil plate	●LDPE (low density polyethylene)

Article code	Article name	Material
	Cardboard 249 x 148 x 100 mm	●Corrugated cardboard, color brown ●Printing FK Logo: Blue- Pantone 300U
	Compound foil 30/70 µm	●PA (Polyamide)/PE (Polyethylene)
	Automat reel without print	●Paper
	Paper 148 x 280 mm	●Paper
Sterilization		
	Indicator tape for EO sterilization	●Process indicator tape for Ethylene Oxide (EO) sterilization
	Sterilization gas	●100% Ethylene oxide